


OPTISPLINT[®] CLINICAL WORKFLOW


SURGICAL

 **PTISPLINT**
by DIGITAL ARCHES



BruxZir[®]
FULL-ARCH IMPLANT PROSTHESIS






Digital Workflow for Full-Arch Prosthesis


Glidewell is proud to partner with Digital Arches to offer OPTISPLINT® to restore full-arch implants **predictably** and **efficiently**.

For BruxZir® Full-Arch Implant Restorations that include ti-bases, **Glidewell recommends the OPTISPLINT workflow.**

GETTING TO KNOW OPTISPLINT®




OPTISPLINT SCAN BODY



FRAME



OPTIWELD™



SCANPLATE™

SCANS NEEDED AT EACH WORKFLOW STAGE

Scans are needed throughout the restorative process.



STAGE 0: PRESURGICAL

- CBCT
 - Dual scan (Fully or partially edentulous cases)
- Intraoral scan (upper/lower)
- Intraoral scan (bite: desired VDO)
- Facial photos or scans



STAGE 1: SURGICAL

Prior to Teeth Extraction

- Intraoral scan of working arch with fiducial marker in place (e.g. ARCHTRACER/teeth/etc.)
- Intraoral scan (opposing arch)
- Intraoral scan (bite: desired VDO)

After Implant Placement

- Tissue scan with fiducial marker (ARCHTRACER™) and Multiunit SCANCAPS
- Implant scan: luted OPTISPLINT on SCANPLATE



STAGE 2: RESTORATIVE


- Facial photos or scans
- Intraoral scan (seated prosthesis)
- Intraoral scan (opposing)
- Intraoral scan (bite: desired VDO)
- Tissue scan with Multiunit SCANCAPS
- Implant scan (luted OPTISPLINT on SCANPLATE)
- Extraoral 360° scan of prosthesis with OPTISPLINT scan analogs

STAGE 0: PRESURGICAL PHASE

The following records are essential for treatment planning full-arch restorations. These records enable the initial design of a prosthesis to align with day-of-surgery scans to fabricate an immediate prosthesis.

Presurgical Records Needed

- CBCT scans
(Dual scan for fully or partially edentulous cases)
- Intraoral scan (upper/lower)
- Intraoral scan (bite: desired VDO)
- Facial photos or scans



Begin case in *My Account*. Indicate day of surgery, type of prosthesis (FP1, FP2, FP3), and upload scans. Glidewell will prepare initial design to be finalized on day of surgery.

STAGE 1: SURGICAL PHASE

THE OUTCOME OF THIS PHASE IS A TEMPORARY PROSTHESIS.

STAGE 1 SCANS NEEDED

The following scans should be acquired on the day of surgery after fiducial marker is placed.

Prior to Tooth Extraction

- Intraoral scan of working arch with fiducial marker in place (e.g. ARCHTRACER/teeth/etc.)
- Intraoral scan (opposing arch)
- Intraoral scan (bite: desired VDO)

After Implant Placement

- Tissue scan: fiducial marker (ARCHTRACER) and Multiunit SCANCAPS
- Implant scan: luted OPTISPLINT on SCANPLATE

SURGICAL WORKFLOW

1 Place Fiducial Markers and Scan

Place fiducial marker (e.g. ARCHTRACER, foundation guide). Take a full-arch scan with intraoral scanner capturing soft tissue and fiducial marker(s).

DENTATE CASES:

Keep preselected teeth as additional markers if teeth are present in the presurgical records.



2 Scan Opposing and Bite

With fiducial marker in place, scan opposing arch and bite.

3 Implant Surgery

Using preferred surgical method, extract desired teeth and place implants. Secure MUAs to each implant.

4 Assemble OPTISPLINT and Lute

Place screws in the OPTISPLINT Scan Body.

With fiducial marker still in place, secure an OPTISPLINT Scan Body on each MUA.





Align the OPTISPLINT wings so that they point toward each other. Lightly hand tighten.

Choose the frame that best reaches all the scan body wings. You may trim the OPTISPLINT Frame, if needed, with provided frame clippers. Rotate the scan bodies as needed to support the frame. Finish hand tightening screws.

Load OPTIWELD(ER)[™] Mixing Gun with the OPTIWELD Dual-Cure Composite cartridge and attach the dispensing tip with the small channel delivery component so that you can flow material onto the OPTISPLINT Frame and down through chimneys.


Working one scan body at a time, flow material through chimney and onto wing and frame. Cure with light, then move to next scan body. Ensure frame is luted securely to all scan bodies.

Make sure you do not get material on the actual scan portion of the sandblasted scan body.



5 Remove OPTISPLINT

Once each scan body has been locked to the frame, unscrew and remove.




6 Tissue Scan

If applicable, remove foundation/bone reduction guide (keep ARCHTRACER in place).

Seat scancaps onto each MUA. Suture and scan using intraoral scanner.


NOTE: Multiunit SCANCAPS utilize the same screws you used with the OPTISPLINT scan bodies. Attach scancaps prior to suturing.



7 Extraoral Scans

Place OPTISPLINT on SCANPLATE and scan.

NOTE: Place OPTISPLINT assembly in enclosed retainer box. Keep in office for Restorative Phase.



8 Send Scans to Glidewell to Fabricate Immediate Provisional



Once Glidewell confirms scans are received, the fiducial markers can be removed.

Glidewell finalizes the design by merging the three surgical day scans with the presurgical planning design to finalize the temporary prosthesis with two options:



Same-Day Delivery: Receive STL file for in-office printing or milling.

Next-Day Delivery: Glidewell overnights temporary (printed or milled) for next-day delivery.

STAGE 2: RESTORATIVE PHASE

THE OUTCOME OF THIS PHASE IS A FIXED TRY-IN PROSTHESIS.

STAGE 2 SCANS NEEDED

The following scans should be acquired in the Restorative Phase. Updated scans are required, even if scans were acquired in Surgical Phase.

- Facial photos or scans
- Intraoral scan (seated prosthesis with surrounding tissue)
- Intraoral scan (opposing arch)
- Intraoral scan (bite: desired VDO)
- Tissue scan with Multiunit SCANCAPS in place
- Extraoral 360° scan of temporary prosthesis with Analog SCANCAPS
- Implant scan: luted OPTISPLINT on SCANPLATE


NOTE: After scans are complete, the physical OPTISPLINT will be sent to Glidewell for final prosthesis fabrication.

RESTORATIVE WORKFLOW

NOTE: Ensure all MUAs are torqued to manufacturer's recommendation prior to seating the OPTISPLINT scan bodies.

1 Scan Provisional Intraorally

Scan seated provisional, ensuring to capture surrounding tissue.



2 Scan Opposing and Bite

Adjust bite if needed.



3 Try in OPTISPLINT

If the fit is not passive, section the OPTISPLINT, secure each section to the MUAs and re-lute.



4 Remove OPTISPLINT


Once each scan body has been locked to the frame, unscrew and remove.



5 Tissue Scan


Seat Multiunit SCANCAPS onto each MUA and scan using intraoral scanner.

NOTE: Multiunit SCANCAPS utilize the same screws you used with the OPTISPLINT scan bodies.




6 Implant Scan (Luted OPTISPLINT on SCANPLATE)

Place OPTISPLINT on SCANPLATE and scan with intraoral scanner.



7 Extraoral 360° Scan of Temporary Prosthesis

Press Analog SCANCAPS into the intaglio side of the temporary prosthesis and scan.



8 Send Files to Lab

Reseat provisional and send files to Glidewell through *My Account*.

FINAL STAGE

1 Send OPTISPLINT to Glidewell

When a try-in is approved, send final scans and OPTISPLINT to Glidewell with the second set of analog scancaps for Glidewell to mill a final prosthesis.



BruxZir®
esthetic
FULL-ARCH IMPLANT PROSTHESIS

Lifetime Warranty and No-Fault Remake
GUARANTEED!

 **Glidewell**

for the sake of smiles