

WHAT'S INCLUDED WITH A MULTI-LEVEL SURGICAL GUIDE?



Implant Verification Prosthesis ■

The implant verification prosthesis (IVP) is connected to the titanium cylinders on the day of surgery and then kept by the restorative dentist. The IVP can be utilized when it is time to fabricate the definitive prosthesis.

Immediate Fixed Screw-Retained Provisional Implant Prosthesis ■ ■

The immediate provisional implant prosthesis is luted to the titanium cylinders on the day of surgery and delivered to the patient to wear until the implants are ready for fabrication of the definitive prosthesis.

Temporary Cylinders ■ ■

The temporary cylinders are connected to the multi-unit abutments and then luted to the immediate provisional implant prosthesis and implant verification prosthesis.

Multi-Unit Abutments ■ ■

Available in straight or angled, 17 or 30 degrees, these abutments seat between the implants and the prosthesis. They allow for the correction of implant angulation to improve screw-access hole locations. They can assist in leveling the restorative platform for easier restoration design and soft-tissue maintenance.

Prosthetic Delivery Guide ■ ■

The prosthetic delivery guide has direction indicators for easy placement of multi-unit abutments. The pegs on the delivery guide allow for the positioning of the immediate provisional prosthesis and the implant verification prosthesis.

Custom Gasket for Relining Prostheses ■

The custom gasket is designed to fit snugly around each of the multi-unit abutments, laying underneath the prosthetic delivery guide. This gasket helps prevent the restorative pick-up material from contaminating the surgical site.

Osteotomy Guide ■ ■ ■

The osteotomy guide allows for the placement of the implants in the exact positions determined by the digital treatment plan. The guide is compatible with any implant system that utilizes a guided drill system.

Mounting Guide ■ ■ ■

The mounting guide allows for the placement of the foundation guide. Pre-designed windows in the mounting guide allow for visualization of the underlying bone when securing the foundation guide in place.

Foundation (Bone Reduction) Guide ■ ■ ■

This is the foundation for the surgery. It is delivered to the bone via a bone-supported mounting guide. It is the foundation upon which all the other guides will be seated. Once secured in place, the flat plane of the foundation guide is utilized to reduce the bone to the predetermined level.

Anchor Pins ■

A 1.5 mm diameter drill is utilized to create the osteotomies for the placement of anchor pins, which are then pushed into position to secure the foundation guide to the bone.

Interlocking Pins ■ ■ ■

These pins secure the multi-level components together. For ease of transition between the various multi-level componentry, the pins are manually placed and released. Floss holes are included on the handles for safety measures.

3D-Printed Segmented Bone Reduction Model ■ ■ ■

This model demonstrates the planned bone reduction.

3D-Printed Model ■ ■ ■

This model demonstrates the bony anatomy of the arch being treated (■ ■ ■) and the post-bone reduction implant placement (■ ■ ■).

■ Tier 1 ■ Tier 2 ■ Tier 3

See reverse for a list of what's included with each tier.

Tier 1	Tier 2	Tier 3
<p>Multi-Level Surgical Guide with Treatment Plan</p>	<p>Multi-Level Surgical Guide with Treatment Plan and Immediate Provisional</p>	<p>Ultimate Multi-Level Surgical Guide Package</p>
<p>What's included:</p> <ul style="list-style-type: none"> ✓ Treatment Plan ✓ 30-Minute Telephone Review Session ✓ 3D-Printed Model ✓ Foundation (Bone Reduction) Guide with Mounting Guide ✓ Osteotomy Guide ✓ Interlocking Pins 	<p>What's included:</p> <ul style="list-style-type: none"> ✓ Treatment Plan ✓ 30-Minute Online Review Session ✓ 3D-Printed Model with Segmented Bone Reduction Model ✓ Foundation (Bone Reduction) Guide with Mounting Guide and Interlocking Pins ✓ Osteotomy Guide ✓ Prosthetic Delivery Guide ✓ Four Multi-Unit Abutments ✓ Four Stock Temporary Cylinders ✓ Immediate Fixed Screw-Retained Provisional Implant Prosthesis 	<p>What's included:</p> <ul style="list-style-type: none"> ✓ Treatment Plan ✓ 30-Minute Online Review Session ✓ 3D-Printed Modular Model with Analogs and Segmented Bone Reduction Model ✓ Foundation (Bone Reduction) Guide with Mounting Guide and Interlocking Pins ✓ Anchor Pins ✓ Osteotomy Guide ✓ Prosthetic Delivery Guide ✓ Multi-Unit Abutments ✓ Custom Gasket for Relining Prostheses ✓ One Set of Stock Temporary Cylinders ✓ Implant Verification Prosthesis ✓ One Set of Custom Precut Temporary Cylinders ✓ Immediate Fixed Screw-Retained Provisional Implant Prosthesis ✓ One Set of Blockout Shims
<p style="text-align: center;">\$1,395*</p> <p style="text-align: center;">Includes 2 Implant Sites</p>	<p style="text-align: center;">\$3,100*</p> <p style="text-align: center;">Includes 4 Implant Sites</p>	<p style="text-align: center;">\$4,950*</p> <p style="text-align: center;">Unlimited Implant Sites</p>

*Price is for implant sites within the same arch and does not include shipping or applicable taxes. Flat-rate price available for BIOMET 3i™ Certain®; CAMLOG® SCREW-LINE; DENTSPLY Implants ANKYLOS® C/X, ASTRA TECH Implant System® and ASTRA TECH Implant System® EV; Hahn™ Tapered Implant System; HIOSSEN® HG System; Inclusive® Tapered Implant System; MegaGen AnyRidge® Implant System; Nobel Biocare Brånemark System® RP, NobelActive® and NobelReplace®, Straumann® Bone Level and Tissue Level; and DentalScrew-Vent®. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components.



Start Your Case Today!

866-497-3692 • glidewell.com/dtp

