

BruxZir[®] **NOW** Pocket Guide

Milling Blocks



BruxZir[®] NOW — Prep and Cement BruxZir Crowns in Just One Appointment!

Specifications

- BruxZir NOW milling blocks are fully sintered zirconia blocks indicated for the chairside production of premolar and molar crowns.
- Sold in packs of three blocks and three burs, BruxZir NOW milling blocks are available exclusively for the TS150[™] Chairside Milling Solution.
- The TS150 mill includes FastDesign[™] software for quick and accurate crown proposals, with finished restorations requiring little to no adjustment after milling.
- Available in 14 shades that correspond to the VITA[®] Classical shade system. (Note: These shades include A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, D2, D3)



Benefits

- **No Oven Wait** of any kind required. Crowns can be designed and milled in less than an hour, with no need for sintering, staining or glazing.
- BruxZir® NOW zirconia provides a strong, tooth-colored alternative to posterior metal occlusal PFMs or cast gold restorations in cases with minimal occlusal space.
- The exceptional marginal integrity and chip-resistant durability of BruxZir NOW zirconia makes it ideal for bruxers who have broken natural teeth or previous restorations.
- Best-in-class chairside features: straightforward preparation, convenient cementation and enduring quality at minimal material thickness.
- With a flexural strength greater than 800 MPa, more than twice that of lithium disilicate, BruxZir NOW is the toughest restoration you can deliver to patients in just one appointment!

BruxZir® NOW

Milling Blocks