

Removable

INCLUSIVE® TRS

LOCATOR BAR OVERDENTURE

Product Highlights

- Implant-retained removable prostheses
- Distribution of occlusal loads reduces implant stress
- Improved functionality with lifelike esthetics
- Predictable implant lab fees and no hidden costs



Description

Our removable implant solutions provide you with reliable restorations that utilize the latest technology for a reasonable and consistent price. Implant overdentures are indicated for patients who desire a more stable prosthesis that improves retention, function, and speech. The **Inclusive Locator® Bar Overdenture** features a CAD/CAM milled bar for distributing occlusal loads more evenly across four or more implants, making it ideal for patients with a heavy bite or when the implants have been placed in softer bone.

Complete Case Includes

- 6 x Inclusive® Tapered Implants



- Final surgical drill



- 6 x Titanium healing abutments



- 6 x Impression copings



- 6 x Analogs



- 4 x Locator® Abutments with Locator® Attachments and processing caps



- Model work

- CAD/CAM milled titanium bar

- Final overdenture with premium denture teeth

Locator is a registered trademark of Zest Anchors.