

Removable

INCLUSIVE® TRS

MINI IMPLANT OVERDENTURE

Product Highlights

- Minimally invasive surgical protocol
- Immediate loading in the presence of primary stability
- Case support from preoperative planning to final prosthesis
- Predictable implant lab fees and no hidden costs



Description

“Mini” or small-diameter implants are traditionally prescribed for patients who do not receive root-form implants due to medical, anatomical or financial reasons. Formerly reserved for transitional use during the osseointegration phase of root-form implants, mini implants have gained acceptance for the immediate loading and long-term stabilization of dentures using a simple and effective O-ring retention system. The **Inclusive Mini Implant Overdenture** provides everything you need for your denture stabilization case: four Inclusive Mini Implants, surgical drill, impression copings, O-ring housings, analogs, model work, and final overdenture — all for a single low price. Additional mini implants, as well as all tools and instruments needed to place them, are also available.

Complete Case Includes

- 4 x Inclusive® Mini Implants
- Cortical bone drill
- 4 x Mini implant impression copings
- 4 x Mini implant analogs
- 4 x O-ring housings (each with additional O-ring)
- Model work
- Final overdenture with Kenson® Teeth



Kenson is a registered trademark of Myerson.