

INCLUSIVE® TRS

SCREW-RETAINED HYBRID DENTURE

Product Highlights

- Precision-milled implant solution
- Fixed-removable restorative option
- All-inclusive pricing with no hidden costs
- Highly functional with lifelike esthetics



Description

The **Inclusive Screw-Retained Hybrid Denture** is a proven replacement for removable prostheses worn by edentulous patients, offering excellent predictability, esthetics, and function. A custom acrylic framework is fabricated on your master model work and then scanned. The scanned data of the framework and the individual implant platforms are transferred to a design station where the one-piece bar framework is finalized. After design, it is milled from titanium for a precise, passive fit. This bar is incorporated into the acrylic processing of the denture tooth setup.

Complete Case Includes

- 6 x Inclusive® Tapered Implants 
- Final surgical drill 
- 6 x Titanium healing abutments 
- 6 x Impression copings 
- 6 x Analogs 
- Model work
- CAD/CAM milled titanium bar
- Final overdenture with premium denture teeth