

Astron CLEARsplint™



	Pearl White		Misty Lavender Glitter
	Pearl White glitter		Misty Lavender
	Brilliant Yellow Glitter		Pumpkin Orange Glitter
	Brilliant Yellow		Pumpkin
	Flamingo Pink Glitter		Red Glitter
	Flamingo Pink		Red
	Green Glitter		Turquoise Blue Glitter
	Green		Turquoise Blue

Nightguards made of the Astron material can have a color tint added, if requested. The appliance will still be translucent, but have an added tint.



**GLIDEWELL
LABORATORIES**

800-407-3326
glidewell.com

.....
Recommends CLEARsplint™

CLEARsplint is a trademark of Astron Dental

GL-2686-011817

CLEARsplint™

THE NIGHTGUARD WITH FLEX!



ADVANTAGES

for both dentist and patient

CLEARsplint is the superior, and optically clear nightguard with **FLEX**

Dentist

- Astron's flexible CLEARsplint is self-adjusting. This superb accuracy offers the closest fit.
- Better fit reduces postoperative chair time
- Dimensional stability means fewer, if any, adjustments
- Composite formulation, a mixture of powder and liquid, can be press-packed or poured
- Ramps and guide planes can be added to reposition the mandible
- Chairside repairs can be made with cold cure acrylics

Patient

- Optical clarity, making it nearly invisible while wearing
- Material is stronger with greater stability, which is ideal for bruxism patients
- Amine-free to never yellow
- MMA*-free to eliminate allergy problems
- Self-adjusting when placed in warm water before wearing
- Wearability and patient comfort are rated excellent

* MMA = methyl methacrylate



CLEARsplint™

CLEARLYdifferent