

• **FIXED RESTORATIVE Rx** •



GLIDEWELL LABORATORIES

4141 MacArthur Blvd. • Newport Beach, CA 92660

800-854-7256 • Fax 800-411-9722
glidewell dental.com




Dr. Name _____ Phone # _____

Acct. # _____ Patient Name _____
First Last



Address/Email _____ **Deliver by 5 p.m. on** See Reverse for Working Times

Enclosed with case: Impressions Models Bite Photos Other: _____

Rx Indicate implant system _____



OCCLUSAL STAINING
 None Light* Medium Dark



**Standard unless specified otherwise*

COMPOSITE RESTORATIONS
 Lava Ultimate Camouflage Composite
 Fiber Reinforcement

IF NO OCCLUSAL CLEARANCE
 Call doctor Spot opposing Metal occlusion
 Would you like this to be a permanent note in your master file?
 Yes No

Signature _____ License # _____
 (see reverse for limited warranty details)

OBSIDIAN PRESSED TO METAL
 Non-Precious* White Noble White High Noble

ZIRCONIA RESTORATIONS
 BruxZir Full-Strength* PrismaTik CZ Lava Plus
 BruxZir Anterior (stump shade recommended for restorations less than 1.5 mm thick) Lava

ALL-CERAMIC RESTORATIONS
 IPS e.max CAD* (Posterior) IPS e.max Press* (Anterior)
 Obsidian All-Ceramic IPS e.max veneer
 Vivaneers No-Prep Veneers*

Indicate stump or present tooth shade for all-ceramics

CUSTOM IMPLANT ABUTMENTS

Choose Material	Choose Brand
<input type="checkbox"/> Titanium*	<input type="checkbox"/> Inclusive*
<input type="checkbox"/> Gold-Colored Titanium	<input type="checkbox"/> BIOMET 3i Encode
<input type="checkbox"/> Zirconia w/ Ti-Base	<input type="checkbox"/> Straumann Variobase
<input type="checkbox"/> Gold Alloy	<input type="checkbox"/> Nobel Biocare NobelProcera

See reverse for supported implant systems

FULL-CAST RESTORATIONS
 Noble-Cast 45 YN (40% Au) OcclusalGold YHN (73.8% Au)
 Noble-Cast 60 YHN (57.5% Au)* JRVY YHN (77% Au)
 Noble-Cast 67 YHN (64% Au) Post & Core

PROVISIONAL RESTORATIONS

Splinted* Individual units
 Transition C&B Smile Transitions

Abutment #(s) _____ Total units _____
 Pontic #(s) _____

BioTemps Provisionals
 Reinforcement: None Wire* Fiber
 BioTemps with Cast-Metal Substructure

Amount of prep reduction: 1 mm* 2 mm
 Perio treatment: Prepare tooth below gingival on tooth #(s) _____ by _____ mm
 Pontic site healing: Prepare ovate socket on tooth #(s) _____ by _____ mm

IN-LAB WORKING TIMES

Please allow full working time for each product selected. Working times are **NOT** guaranteed and do **NOT** include weekends or holidays. Rush service available on most products but must be prescheduled (see below).

IN-LAB WORKING TIMES


Obsidian Pressed to Metal	5 days	Partials & Dentures	
Diagnostic Wax-Up	5 days	Metal Frames	6 days
Obsidian w/ attachments	7 days	Metal Frame with teeth and wax.....	11 days
Custom Abutments	8 days	Metal Partial to finish	11 days
All-ceramic/Zirconia restorations	5 days	Occlusion rims/Custom tray	3 days
Full-cast restorations	5 days	Denture setup try-in	5 days
Composite restorations	5 days	Denture try-in to finish	5 days
Camouflage Composite.....	3 days	Valplast/tcs setup teeth in wax.....	5 days
Fiber-reinforced composites	5 days	Valplast/tcs start to completion.....	7 days
BioTemps Provisionals	5 days	Denture soft liner/acrylic repair or reline.....	3 days
With cast-metal substructure	6 days	Partial repair/reline.....	5 days
With cast-metal substructure over implant.....	6 days	PlaySafe mouthguards	4 days
Screw-retained over implant	6 days	Silent Nite sl	4 days
Smile Transitions appliance	7 days	TAP/TAP 3/EMA	5 days
Transition Crowns & Bridges	6 days	Comfort H/S Bite Splint	4 days
Lava	5 days	Nightguards/Bleach trays/Custom trays	4 days

All rush cases must be prescheduled


by calling **800-944-7874** before the case is shipped. Time of pickup and delivering may affect turnaround time.

PREPARATION GUIDELINES


OBSIDIAN PFM ANTERIOR




OBSIDIAN PFM POSTERIOR




OBSIDIAN PFM-PORCELAIN LABIAL OR 360° MARGIN




ALL-CERAMIC/COMPOSITE VENEERS




INLAY



ONLAY



ALL-CERAMIC/COMPOSITE CROWNS



TERMS AND WARRANTY INFORMATION

All restorations made in the USA.
Only \$7 shipping per box EACH way (contiguous U.S. only; shipping charge varies for Alaska, Hawaii and Puerto Rico).

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

LIMITED WARRANTY/LIMITATION OF LIABILITY. Glidewell Laboratories ("the lab") warrants that all dental devices (a "device") are made according to your specification and approval in the belief that the device will be useful and MAKES NO OTHER WARRANTIES INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Subject to the return of a device that is placed and then fails, the lab will repair or replace the device without charge for the cost of materials and workmanship or refund the original price paid, at the lab's option, as follows: (1) screw-retained implant abutments up to 20 years (including replacement of a device that was made by the lab that was attached to said abutments and including an Inclusive® Tapered titanium dental implant upon request at no extra charge); (2) porcelain to metal, all ceramic, all metal, single-unit inlay, onlay and crown composite resin final prosthetics, and screw-retained crowns and implant bars, up to seven years; (3) composite resin bridges (excluding Maryland and inlay/onlay bridges) up to five years; (4) Transition Crowns and Bridges® up to two years; (5) dentures and including screw-retained dentures but excluding immediate and provisional dentures up to one year if the failure is due to defects in materials or workmanship; (6) thermoformed appliances and splints if the failure is due to defects in materials or workmanship, composite resin, Maryland and inlay/onlay bridges, flippers and Reflex™ Plus and provisional partials up to six months; (7) Smile Transitions™ cosmetic appliances, immediate and provisional dentures, partials, retainers, surgical stents and radiographic guides, and all other dental devices up to thirty days if the failure is due to defects in materials or workmanship. This Warranty is exclusively for your benefit, is not transferable and does not extend to any patients. You agree to pay all other costs of adjustment, repair and replacement of a device. Except where prohibited by law, the lab WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES ARISING FROM THE USE OF A DEVICE, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL, regardless of the theory asserted, including warranty, contract, negligence or strict liability and if such disclaimer is not permitted by law, the duration of any implied warranty is limited to 90 days from the date of delivery. In the event of a dispute and absent an amicable resolution the parties mutually agree to waive class actions in favor of mandatory individual arbitration of claims under this limited warranty in and in accordance with the laws of California. The lab does not guarantee the performance of independent carriers.

INCLUSIVE® TITANIUM AND ZIRCONIA W/TI-BASE CUSTOM ABUTMENTS ARE COMPATIBLE WITH THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™ Certain®	CAMLOG® SCREW-LINE	DENTSPLY Implants ANKYLOS® C/X ASTRA TECH Implant System®	Glidewell Direct Hahn™ Tapered Implant System Inclusive® Tapered Implant System
HIOSSSEN® HG System	Keystone Dental PrimaConnex®	Neoss® Neoss®	Nobel Biocare Brånemark System® RP NobelActive® NobelReplace®
Straumann® Bone Level Tissue Level	Sweden & Martina Premium Shelta	Zimmer Dental Screw-Vent®	

Inclusive is a registered trademark of Prismatic Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prismatic Dentalcraft, Inc. All other trademarks are property of their respective owners. Glidewell Laboratories works in partnership with Keystone Dental, Neoss and Sweden & Martina.