

• **BioTemps Provisionals Rx** •



**GLIDEWELL  
LABORATORIES**

4141 MacArthur Blvd. • Newport Beach, CA 92660

**800-777-0695 • Fax 800-411-9722**  
**glidewell.com**



Dr. Name \_\_\_\_\_ Phone # \_\_\_\_\_

Acct. # \_\_\_\_\_ Patient Name \_\_\_\_\_  
First Last

Address/Email \_\_\_\_\_ **Deliver by 5 p.m. on** See Reverse for Working Times

Enclosed with case:  Impressions  Models  Bite  Photos  Other: \_\_\_\_\_

Rx




Indicate Shade Here

**SELECT INCISAL EDGE TYPE**



Rounded  Squared  Pointed

**SELECT PONTIC DESIGN TYPE**



Pontic site healing: Prepare ovate socket on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_mm

Splinted\*  Individual units  
 Cement-on implant  Screw-retained implant

Abutment #(s) \_\_\_\_\_

Pontic #(s) \_\_\_\_\_

Extract #(s) \_\_\_\_\_

Extend Gingival Margins:  
Tooth #(s) \_\_\_\_\_, \_\_\_\_\_mm

Amount of Prep Reduction:  1 mm\*  2 mm

Perio treatment: Prepare tooth below gingival on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_mm

Perio site healing: Prepare ovate socket on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_mm

**CHOOSE RESTORATION**

**BioTemps Provisionals**  
• 6-month warranty / 5 working days  
Reinforcement:  None  Wire\*  Fiber

**BioTemps with Cast-Metal Substructure**

**Transition Crowns & Bridges**

*(No reinforcement available)*

*Transition C&B are milled from an acetal resin block*

- 2-year warranty / 6 working days
- Fabricated using a virtual CAD preparation
- Design files are stored for 2 years

Fabricate Over:

Lab Prep  Solid Model  
 Die-Trim Model  Implant

**Smile Transitions**

*(No reinforcement available)*

*Smile Transitions are milled from an acetal resin block*

- 30-day warranty / 7 working days
- Fabricated as a no-prep temporary solution
- Design files are stored for 2 years

Upper:  
Tooth # \_\_\_\_\_ to Tooth # \_\_\_\_\_

Lower:  
Tooth # \_\_\_\_\_ to Tooth # \_\_\_\_\_

Occlusion:

Maintain Occlusion\*  
 Increase VDO \_\_\_\_\_mm

*\*Standard unless specified otherwise*

Signature \_\_\_\_\_ License # \_\_\_\_\_ Date \_\_\_\_\_

(see reverse for limited warranty details)

**IN-LAB WORKING TIMES**

All restorations made in the USA.

Please allow full working time for each product selected. Working times are **NOT** guaranteed and do **NOT** include weekends or holidays. Rush service available on most products but must be prescheduled (see below).


**IN-LAB WORKING TIMES**

<b>Obsidian/Diagnostic wax-up</b> .....	5 days	<b>Partials &amp; Dentures</b>	
<b>Obsidian w/ attachments</b> .....	7 days	Metal Frames .....	6 days
<b>Custom abutments</b> .....	8 days	Metal Frame with teeth and wax.....	11 days
<b>All-ceramic/Zirconia restorations</b> .....	5 days	Metal Partial to finish .....	12 days
<b>Full-cast restorations</b> .....	5 days	Occlusion rims/Custom tray .....	3 days
<b>Composite restorations</b> .....	3 days	Denture setup try-in .....	5 days
<b>Fiber-reinforced composites</b> .....	5 days	Denture try-in to finish .....	5 days
<b>BioTemps Provisionals</b> .....	5 days	Valplast/tcs setup teeth in wax.....	5 days
With cast-metal substructure .....	6 days	Valplast/tcs start to completion .....	7 days
With cast-metal substructure over implant ..	6 days	Denture soft liner/acrylic repair or reline.....	5 days
Screw-retained over implant.....	6 days	Partial repair/reline .....	5 days
<b>Smile Transitions appliance</b> .....	7 days	<b>PlaySafe mouthguards/Silent Nite sl</b> .....	4 days
<b>Transition Crowns &amp; Bridges</b> .....	6 days	<b>TAP/TAP 3/EMA</b> .....	5 days
<b>Lava</b> .....	5 days	<b>Comfort H/S Bite Splint</b> .....	4 days
		<b>Nightguards/Bleach trays/Retainers</b> .....	4 days

All rush cases must be prescheduled by calling **800-944-7874** before the case is shipped. Time of pickup and delivery may affect turnaround time.


**PREPARATION GUIDELINES**

**OBSIDIAN TO METAL ANTERIOR**




- A. 2 mm incisal reduction
- B. 1.5 mm middle third reduction
- C. Labial and lingual walls must be convergent.
- D. Preparation should be cut in three planes.

**OBSIDIAN TO METAL POSTERIOR**




- A. 2 mm occlusal reduction
- B. 1.5 mm middle third reduction
- C. Buccal and lingual walls must be convergent
- D. Preparation should be cut in three planes

**OBSIDIAN TO METAL – CERAMIC LABIAL OR 360° MARGIN**




1.25 mm gingival reduction using rounded shoulder margin design

**ALL-CERAMIC/COMPOSITE VENEERS**




- A. 0.3 to 1.0 mm labial reduction

**INLAY**




- A. 1.5 to 2 mm occlusal reduction
- B. Round all sharp line angles and occlusal edges, and eliminate undercuts.
- C. Proximal and occlusal walls should have 6-8 degrees taper

**ONLAY**




**ALL-CERAMIC/COMPOSITE CROWNS**



**TERMS AND WARRANTY INFORMATION**

Only \$7 shipping per box each way (contiguous U.S. only; shipping charge varies for Alaska, Hawaii and Puerto Rico).



**All Restorations Made in the USA**

***We honor VISA, MASTERCARD, AMEX and DISCOVER.***

**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** Glidewell Laboratories is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [glidewell.com/policies-and-warranties/](http://glidewell.com/policies-and-warranties/).

**FOR LAB USE ONLY**  
TELEPHONE CALL RECORD

DR. \_\_\_\_\_ ACCT. # \_\_\_\_\_

PATIENT NAME \_\_\_\_\_

RE: \_\_\_\_\_

RESULT \_\_\_\_\_

**DATE DUE IN OFFICE** \_\_\_\_\_

DATE OF CALL \_\_\_\_\_ INITIALS \_\_\_\_\_