

• **FIXED RESTORATIVE Rx** •



**GLIDEWELL  
LABORATORIES**

4141 MacArthur Blvd. • Newport Beach, CA 92660

**800-854-7256 • Fax 800-411-9722**  
**glidewell dental.com**


Dr. Name \_\_\_\_\_ Phone # \_\_\_\_\_

Acct. # \_\_\_\_\_ Patient Name \_\_\_\_\_  
First Last

Address/Email \_\_\_\_\_ **Deliver by 5 p.m. on** See Reverse for Working Times

Enclosed with case:  Impressions  Models  Bite  Photos  Other: \_\_\_\_\_

**Rx** Implant System \_\_\_\_\_ Implant Diameter \_\_\_\_\_ mm



Signature \_\_\_\_\_ License # \_\_\_\_\_  
(see reverse for limited warranty details)

**FINAL CERAMIC SHADE**



Indicate Shade Here \_\_\_\_\_

**PRESENT TOOTH OR STUMP SHADE**



Indicate Shade Here \_\_\_\_\_


**OCCLUSAL STAINING**

None  Light\*  Medium  Dark

**PONTIC DESIGN**



**MARGIN AND METAL DESIGN**



\*Standard unless specified otherwise

**COMPOSITE RESTORATIONS**

Camouflage NanoHybrid Composite  
 Fiber Reinforcement

**IF NO OCCLUSAL CLEARANCE**

Call doctor  Would you like this to be a permanent note in your master file?  
 Spot opposing  Yes  No  
 Metal occlusion  Yes  No

**OBSIDIAN FUSED TO METAL**

Non-Precious\*  White Noble  White High Noble

**ZIRCONIA RESTORATIONS**

BruxZir Full-Strength\* (1,150 MPa)  
 BruxZir Anterior (650 MPa) (stump shade recommended for restorations less than 1.5 mm thick)  
 Lava Crowns & Bridges  PrismaTik Clinical Zirconia

**ALL-CERAMIC RESTORATIONS**

IPS e.max CAD\* (Posterior)  IPS e.max Press\* (Anterior)  
 Obsidian All-Ceramic  IPS e.max veneer  
 Vivaneers No-Prep Veneers\*

Indicate stump or present tooth shade for all-ceramics

**CUSTOM IMPLANT ABUTMENTS**

Choose Material Choose Brand

Titanium\*  Inclusive\*  
 Gold-Colored Titanium  BIOMET 3i Encode  
 Zirconia w/ Ti-Base  Straumann Variobase  
 Gold Alloy

See reverse for supported implant systems

**SCREW-RETAINED IMPLANT CROWNS**

BruxZir Full-Strength  BruxZir Anterior  IPS e.max

Specify implant brand, system and diameter on Rx

**FULL-CAST RESTORATIONS**

Noble-Cast 42 YN (2% Au)  OcclusalGold YHN (73.8% Au)  
 Noble-Cast 45 YN (40% Au)  JRVT YHN (77% Au)  
 Noble-Cast 60 YHN (57.5% Au)\*  Post & Core  
 Noble-Cast 67 YHN (64% Au)

**PROVISIONAL RESTORATIONS**

Splinted\*  Individual units  
 Transition C&B  Smile Transitions

Abutment #(s) \_\_\_\_\_ Total units \_\_\_\_\_

Pontic #(s) \_\_\_\_\_ Total units \_\_\_\_\_

BioTemps Provisionals

Reinforcement:  None  Wire\*  Fiber

BioTemps with Cast-Metal Substructure

Amount of prep reduction:  1 mm\*  2 mm

Perio treatment: Prepare tooth below gingival on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_ mm

Pontic site healing: Prepare ovate socket on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_ mm

**IN-LAB WORKING TIMES**

Please allow full working time for each product selected. Working times are **NOT** guaranteed and do **NOT** include weekends or holidays. Rush service available on most products but must be prescheduled (see below).

**IN-LAB WORKING TIMES**


<b>Obsidian Fused to Metal</b> .....	5 days	<b>Partials &amp; Dentures</b>	
<b>Diagnostic wax-up</b> .....	5 days	Metal Frames .....	6 days
<b>Obsidian w/ attachments</b> .....	7 days	Metal Frame with teeth and wax.....	11 days
<b>Custom abutments</b> .....	8 days	Metal Partial to finish .....	12 days
<b>All-ceramic/zirconia restorations</b> .....	5 days	Occlusion rims/Custom tray .....	3 days
<b>Full-cast restorations</b> .....	5 days	Denture setup try-in .....	5 days
<b>Composite restorations</b> .....	3 days	Denture try-in to finish .....	5 days
<b>Fiber-reinforced composites</b> .....	5 days	Valplast/tcs setup teeth in wax.....	5 days
<b>BioTemps Provisionals</b> .....	5 days	Valplast/tcs start to completion.....	7 days
With cast-metal substructure .....	6 days	Denture soft liner/acrylic repair or reline.....	5 days
With cast-metal substructure over implant .....	6 days	Partial repair/reline.....	5 days
Screw-retained over implant .....	6 days	<b>PlaySafe mouthguards</b> .....	4 days
<b>Smile Transitions appliance</b> .....	7 days	<b>Silent Nite sl</b> .....	4 days
<b>Transition Crowns &amp; Bridges</b> .....	6 days	<b>TAP/TAP 3/dreamTAP/EMA</b> .....	5 days
<b>Lava</b> .....	5 days	<b>Comfort H/S Bite Splint</b> .....	4 days
		<b>Nightguards/Bleach trays/Custom trays</b> .....	4 days

**All rush cases must be prescheduled**

by calling **800-944-7874** before the case is shipped. Time of pickup and delivering may affect turnaround time.


**PREPARATION GUIDELINES**

**OBSIDIAN TO METAL ANTERIOR**




- A. 2 mm incisal reduction
- B. 1.5 mm middle third reduction
- C. Labial and lingual walls must be convergent
- D. Preparation should be cut in three planes

**OBSIDIAN TO METAL POSTERIOR**




- A. 2 mm occlusal reduction
- B. 1.5 mm middle third reduction
- C. Buccal and lingual walls must be convergent
- D. Preparation should be cut in three planes

**OBSIDIAN TO METAL - CERAMIC LABIAL OR 360° MARGIN**




**ALL-CERAMIC/COMPOSITE VENEERS**




- A. 0.3 to 1 mm labial reduction

**INLAY**




- A. 1.5 to 2 mm occlusal reduction
- B. Round all sharp line angles and occlusal edges, and eliminate undercuts.
- C. Proximal and occlusal walls should have 6-8 degrees taper.

**ONLAY**




**ALL-CERAMIC/COMPOSITE CROWNS**



**TERMS AND WARRANTY INFORMATION**

Only \$7 shipping per box EACH way (contiguous U.S. only; shipping charge varies for Alaska, Hawaii and Puerto Rico).



**All Restorations Made in the USA**

*We honor VISA, MASTERCARD, AMEX and DISCOVER.*

**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** Glidewell Laboratories is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [glidewell.com/policies-and-warranties/](http://glidewell.com/policies-and-warranties/).

**INCLUSIVE® CUSTOM TITANIUM AND ZIRCONIA W/TI-BASE ABUTMENTS ARE COMPATIBLE WITH THE FOLLOWING IMPLANT SYSTEMS**

<b>BIOMET 3i™</b> Certain®	<b>CAMLOG®</b> SCREW-LINE	<b>DENTSPLY Implants</b> ANKYLOS® C/X ASTRA TECH Implant System®	<b>Glidewell Direct</b> Hahn™ Tapered Implant System Inclusive® Tapered Implant System
-------------------------------	------------------------------	---	--

<b>HIOSSSEN®</b> HG System	<b>Keystone Dental</b> PrimaConnex®	<b>MegaGen</b> AnyRidge® Implant System	<b>Neoss®</b> Neoss® Implant System
-------------------------------	--	---	---

<b>Nobel Biocare</b> Brånemark System® RP NobelActive® NobelReplace®	<b>Straumann®</b> Bone Level Tissue Level	<b>Sweden &amp; Martina</b> Premium Shelta	<b>Zimmer Dental</b> Screw-Vent®
---	---	--	-------------------------------------

Inclusive is a registered trademark of Prismatic Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prismatic Dentalcraft, Inc. All other trademarks are property of their respective owners. Glidewell Laboratories works in partnership with Keystone Dental, Neoss and Sweden & Martina.