



# GLIDEWELL LABORATORIES

4141 MacArthur Blvd. • Newport Beach, CA 92660

**800-854-7256 • Fax 800-411-9722**  
**glidewell dental.com**

Dr. Name \_\_\_\_\_ Phone # \_\_\_\_\_

Acct. # \_\_\_\_\_ Patient Name \_\_\_\_\_  
First Last

Address/Email \_\_\_\_\_ **Deliver by 5 p.m. on** See Reverse for Working Times

Enclosed with Case:  Impressions  Models  Bite  Photos  Other: \_\_\_\_\_

### PROVISIONAL RESTORATIONS

- BioTemps Provisionals  
 Reinforcement:  None  Wire\*  Fiber  
 Cast-Metal
- Transition C&B  Smile Transitions
- Abutment #(s) \_\_\_\_\_
- Pontic #(s) \_\_\_\_\_ Total units \_\_\_\_\_
- Splinted\*  Cement-on Implant
- Individual Units  Screw-Retained Implant
- Amount of prep reduction:  1 mm\*  2 mm
- Perio treatment: Prepare tooth below gingival on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_ mm
- Pontic site healing: Prepare ovate socket on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_ mm

### DENTURES/FLIPPERS/FLEXIBLE PARTIALS

- Denture  Flipper  Valplast  tcs  DuraFlex
- Custom tray  Bite rim  Wax setup try-in  Finish
- Kenson Teeth (Included at no extra charge)
- Premium Brand Teeth (Extra charge applies)**  
 Shade \_\_\_\_\_ Brand \_\_\_\_\_ Mould \_\_\_\_\_

### NIGHTGUARDS/RETAINERS

- Upper  Lower  Scan/Save File
- Comfort H/S\* (hard/soft)  Astron CLEARsplint (hard processed acrylic)
- Soft nightguard  Comfort (hard)
- Clear-Lock Retainers for Life:  
 set of 2  set of 4  set of 6
- Clear ortho  Hawley  QCM Clear-Wire

### SNORING/SLEEP APNEA APPLIANCES

- Silent Nite sl\*  EMA  Scan/Save File
- Buy 1  dreamTAP
- Buy 2 and save  TAP 3  TAP

### PLAYSAFE MOUTHGUARDS

- Jr  Lt  Lt Pro  Med\*  Hvy  Hvy Pro
- Helmet strap **Specify color(s) on Rx**

### ZIRCONIA/ALL-CERAMIC RESTORATIONS

- BruxZir Full-Strength\* (1,150 MPa)  BruxZir Esthetic (870 MPa)
- Obsidian All-Ceramic  IPS e.max (anterior)
- PrismaTik Clinical Zirconia  IPS e.max (posterior)
- Lava Zirconia

### VIVANEER VENEER

- Conventional-Prep:  
 BruxZir Esthetic  Obsidian  IPS e.max
- Minimal-Prep:  
 Obsidian  IPS e.max

### OBSIDIAN FUSED TO METAL

- Non-Precious\*  White Noble  White High Noble

### SCREW-RETAINED RESTORATIONS

- BruxZir Full-Strength\*  Obsidian Fused to White Noble (w/ Ti-Base)
- BruxZir Esthetic  Obsidian Fused to White High Noble (w/ Ti-Base)
- IPS e.max (w/ Ti-Base)  PrismaTik CZ (w/ Ti-Base)

### CUSTOM ABUTMENTS

- Titanium\*  Gold-Tone Titanium  Gold Alloy
- Zirconia w/ Ti-Base  BIOMET 3i Encode
- Prepare existing abutment

### COMPOSITE RESTORATIONS

- Camouflage NanoHybrid Composite  Fiber Reinforcement

### FULL-CAST RESTORATIONS

- Noble-Cast 42 YN (2% Au)  JRVT YHN (77% Au)
- Noble-Cast 45 YN (40% Au)  Non-Precious
- Noble-Cast 60 YHN (57% Au)\*  White Noble
- Noble-Cast 67 YHN (64% Au)  WHN (45% Au)
- OcclusalGold YHN (73% Au)  Post & Core

### METAL PARTIALS

- Vitallium 2000\*  Vitallium 2000 Plus  tcs/Vitallium
- Valplast/Vitallium  Wironium  tcs/Wironium
- Metal frame try-in  Printed frame try-in  Frame w/occlus. rim
- Frame w/setup try-in  Finish  Scan/Save File

### FINAL SHADE



### STUMP SHADE

*Must indicate prepped tooth shade for all-ceramics*

### OCCLUSAL STAINING

- None  Light\*
- Medium  Dark

### PONTIC DESIGN



### MARGIN AND METAL DESIGN




Labial Butt Junction Junct. Collar Excl. Cusp Incl. Cusp

### IF NO OCCLUSAL CLEARANCE

- Call doctor  Spot opposing  Metal occlusion
- Metal island  Make this a permanent note

**Rx**



Signature \_\_\_\_\_ License # \_\_\_\_\_  
 (see reverse for limited warranty details)

## IN-LAB WORKING TIMES


Please allow full working time for each product selected. Working times are NOT guaranteed and do NOT include weekends or holidays. Rush service available on most products but must be prescheduled (see below).

<b>Obsidian Fused to Metal</b> .....5 days	<b>BioTemps Provisionals</b> .....5 days	Metal Partial to finish ..... 12 days
<b>Diagnostic wax-up</b> .....5 days	With cast-metal substructure .....6 days	Occlusion rims/Custom tray ..... 3 days
<b>Obsidian w/ attachments</b> .....7 days	With cast-metal substructure over implant .....6 days	Denture setup try-in ..... 5 days
<b>Custom abutments</b> .....8 days	Screw-retained over implant .....6 days	Finish after setup try-in..... 5 days
<b>All-ceramic screw-retained restorations</b> .....8 days	<b>Smile Transitions appliance</b> .....7 days	Denture try-in to finish ..... 5 days
<b>Obsidian Fused to Metal screw-retained restorations</b> .....8 days	<b>Transition Crowns &amp; Bridges</b> .....6 days	Valplast/tcs setup teeth in wax ..... 5 days
<b>All-ceramic/zirconia restorations</b> .....5 days	<b>Partials &amp; Dentures</b>	Finish after setup try-in..... 6 days
<b>Full-cast restorations</b> .....5 days	Metal Frames .....6 days	Valplast/tcs start to completion..... 7 days
<b>Composite restorations</b> .....3 days	Metal Frame with teeth and wax .....11 days	Denture soft liner/acrylic repair or reline..... 5 days
<b>Fiber-reinforced composites</b> .....5 days	Finish after setup try-in.....6 days	Partial repair/reline..... 5 days
		<b>PlaySafe mouthguards/Silent Nite sl</b> ..... 4 days
		<b>TAP/TAP 3/dreamTAP/EMA</b> ..... 5 days
		<b>Comfort H/S Bite Splint</b> ..... 4 days
		<b>Nightguards/Bleaching trays/Retainers</b> ..... 4 days

**All rush cases must be prescheduled by calling 800-944-7874 before the case is shipped. Time of pickup and delivery may affect turnaround time.**


## PREPARATION GUIDELINES

**Anterior: BruxZir Full-Strength**




- A. 1 mm ideal reduction (0.6 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Labial and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Anterior: BruxZir Esthetic**




- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or shoulder margins preferred
- C. Labial and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Anterior: Obsidian and IPS e.max**




- A. 1.5 mm ideal reduction (1.25 mm minimum)
- B. Shoulder margins preferred
- C. Labial and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Anterior: Obsidian Fused to Metal**




- A. 2 mm ideal reduction (1.5 mm minimum)
- B. Chamfer or Shoulder margins preferred
- C. Labial and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Posterior: BruxZir Full-Strength**




- A. 1 mm ideal reduction (0.6 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Buccal and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Posterior: Obsidian and IPS e.max**



- A. 1.5 mm ideal reduction (1.25 mm minimum)
- B. Shoulder margins preferred
- C. Buccal and lingual walls must be convergent
- D. Preparation should be cut in three planes

**Posterior: Obsidian Fused to Metal**



- A. 2 mm ideal reduction (1.5 mm minimum)
- B. Chamfer or Shoulder margins preferred
- C. Buccal and lingual walls must be convergent
- D. Preparation should be cut in three planes

## TERMS AND WARRANTY INFORMATION

We honor VISA, MASTERCARD, AMEX and DISCOVER.

**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** Glidewell Laboratories is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [glidewell.com/policies-and-warranties/](http://glidewell.com/policies-and-warranties/).



**All Restorations Made in the USA**