

Dr. Name \_\_\_\_\_ Phone # \_\_\_\_\_

Acct. # \_\_\_\_\_ Patient ID/Name \_\_\_\_\_ First Last

Address/Email \_\_\_\_\_ Deliver by 5 p.m. on \_\_\_\_\_

Enclosed with Case:  Impressions  Models  Bite  Photos  Other: \_\_\_\_\_

PROVISIONAL RESTORATIONS
 BioTemps Provisionals
Reinforcement:  None  Wire  Fiber
 Cast-Metal
 Transition C&B  Smile Transitions
Abutment #(s) \_\_\_\_\_
Pontic #(s) \_\_\_\_\_ Total units \_\_\_\_\_
 Splinted  Cement-On Implant
 Individual Units  Screw-Retained Implant
Amount of prep reduction:  1 mm  2 mm
 Perio treatment: Prepare tooth below gingival
on tooth #(s) \_\_\_\_\_ by \_\_\_\_\_ mm
 Pontic site healing: Prepare ovate socket on tooth
#(s) \_\_\_\_\_ by \_\_\_\_\_ mm

DENTURES/FLIPPERS/FLEXIBLE PARTIALS
Denture Partial Select Phase
 Handcrafted Denture  Flipper  Custom tray
 Digital Denture  Valplast  Bite rim
 Immediate Digital Denture  tcs  Setup try-in
 Immediate Denture  DuraFlex  Finish
 Duplicate Digital Denture
 Reference Digital Denture

Digital Denture Teeth
Shade \_\_\_\_\_ Mould \_\_\_\_\_
Select Teeth for Partials and Handcrafted Dentures
 Kenson Teeth (Standard)
Shade \_\_\_\_\_ Mould \_\_\_\_\_
 Premium Brand Teeth (Extra Charge)
Shade \_\_\_\_\_ Mould \_\_\_\_\_ Brand \_\_\_\_\_

SIMPLY NATURAL METAL PARTIALS
Metal frame with acrylic and Kenson teeth
Frame Material Phase
 SLM-printed cobalt chrome frame  Metal frame try-in
 Printed frame try-in
 Frame w/occlus. rim
 Frame w/setup try-in
 Finish
 Scan/Save File (extra charge applies)
 Lab select complete design

SNORING/SLEEP APNEA APPLIANCES
 Buy 1  Buy 2 and save  Scan/Save File
 Silent Nite 3D (digital impressions only)  Silent Nite
 Silent Nite with Glidewell Hinge  OASYS Hinge Appliance
 EMA  flexTAP  dreamTAP  TAP 3 TL

ZIRCONIA & ALL-CERAMIC RESTORATIONS
 BruxZir Full-Strength™ (>1,000 MPa)  IPS e.max
 BruxZir Radiant (778 MPa)  Bilayered Clinical Zirconia
 BruxZir Esthetic™ (870 MPa)
 BruxZir Fusion (Veneer 778 MPa, Core > 1,000 MPa)

VENEERS
 BruxZir Esthetic veneer\*  IPS e.max veneer  Layered IPS e.max veneer

PFM
 Obsidian to Non-Precious  Obsidian to White High Noble
 Obsidian to White Noble  Porcelain fused to OcclusalGold YHN

FULL-CAST RESTORATIONS
 Noble-Cast 45 YN (40% Au)  White Noble
 Noble-Cast 60 YHN (57.5% Au)  White High Noble (40% Au)
 Non-Precious  Post & Core

SCREW-RETAINED RESTORATIONS
 BruxZir Full-Strength™  Obsidian to White Noble
 BruxZir Esthetic  Obsidian to White High Noble
 IPS e.max  Bilayered Clinical Zirconia
 BruxZir Fusion

SCREWMOUNTABLE
(crown with screw-access hole cemented over custom abutment)
 BruxZir Full-Strength™  BruxZir Esthetic
 BruxZir Fusion  IPS e.max

CUSTOM ABUTMENTS
 Titanium  Gold-Tone Titanium  Zirconia w/ Ti-Base
 BIOMET 3i Encode  Gold Alloy  Prepare existing abutment

FULL-ARCH IMPLANTS
 BruxZir Full-Arch Implant Prosthesis  Screw-Retained Hybrid Denture
 Full-Strength  Esthetic  Locator Overdenture

NIGHTGUARDS/RETAINERS/MIGRAINE PREVENTION
 Upper  Lower
 Buy 1  Buy 2 and save  Scan/Save File
 Comfort3D  Comfort H/S  Comfort Stabilizer (3D-printed, hard) (clear, hard with soft reline) (buy 2 savings not available)
 Soft nightguard (buy 2 savings not available)  Bleaching Tray (buy 2 savings not available)
 CLEARsplint (self-adjusting, hard)  NTL-tss Plus  NTL OmniSplint
 Clear-Lock Retainers:  Qty of 2  Qty of 4  Qty of 6
 Clear ortho  Hawley  Essix Retainer (1 tooth)

PLAYSAFE MOUTHGUARDS
 Jr  Lt  Lt Pro  Med\*  Hvy  Hvy Pro
 Helmet strap \_\_\_\_\_ Specify color(s) on Rx
Name: \_\_\_\_\_

ORB SPORT MOUTHGUARD
 ORB Sport Mouthguard Includes digital file storage of the patient's model for 7 years


IF NO OCCLUSAL CLEARANCE
 Call doctor  Spot opposing  Metal occlusion  Metal island  Make this a permanent note

STUMP SHADE FINAL SHADE
OCCLUSAL STAINING
 None  Light\*  Medium  Dark
PONTIC DESIGN

FOR BRUXZIR FUSION CASES
INCISAL LOBE DESIGN:  Less  Light\*  Heavy  None
INCISAL TRANSLUCENCY:  Less  Light\*  Heavy
ANTERIOR DESIGN STYLE:  Triangle  Round  Square

ANATOMICAL SURFACE TEXTURE
 None  Light\*  Medium

Rx Implant System \_\_\_\_\_ (If applicable)
Implant Diameter \_\_\_\_\_ mm
(See reverse for flat-rate pricing info)



Signature \_\_\_\_\_
License \_\_\_\_\_ Date \_\_\_\_\_

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

\*Standard unless specified otherwise. \*\*BruxZir Esthetic is the standard for anterior restorations. BruxZir Full-Strength is the standard for posterior restorations.

## TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

*We honor VISA, MASTERCARD, AMEX and DISCOVER.*

**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. *Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.* Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [glidewell.com/policies-and-warranties](http://glidewell.com/policies-and-warranties).



• BruxZir Restorations



• Custom Abutments



- All-Ceramic Restorations
- PFM Restorations
- Full-Cast Restorations



- BioTemps® Provisionals
- Bite Splints
- Mouthguards and Nightguards




- Silent Nite
- Silent Nite 3D
- Silent Nite with Glidewell Hinge
- flexTAP



- dreamTAP
- TAP 3 TL
- EMA

*All rush cases must be prescheduled by calling 800-944-7874 before the case is shipped. Time of pickup and delivery may affect turnaround time.*

## PREPARATION GUIDELINES



### BruxZir Radiant and BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

### BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

## FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BioHorizons®  
Tapered Internal

BIOMET 3i™  
Certain®

CAMLOG®  
SCREW-LINE

DENTSPLY Implants  
ANKYLOS® C/X  
ASTRA TECH Implant System®  
ASTRA TECH Implant System® EV

Nobel Biocare  
Brånemark System® RP  
NobelActive®  
NobelReplace®

Glidewell Direct  
Glidewell HT™ Implant System  
Hahn™ Tapered Implant System  
Inclusive® Tapered Implant System

HIOSSEN®  
HG System

Implant Direct®  
InterActive®  
Legacy®  
RePlant®

Keystone® Dental  
PrimaConnex®

Straumann®  
BLX  
Bone Level  
Bone Level SC  
Tissue Level

MegaGen  
AnyRidge® Implant System

MIS®  
C1 Implant System  
SEVEN®

NeoDent®  
GM

Neoss®†  
ProActive® Narrow  
ProActive® Standard

Zimmer Dental  
Eztetic®  
Screw-Vent®

Restorations for the implant systems above are fabricated using components manufactured by PrismaTik Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. †For Neoss®, flat-rate pricing is only available for screw-retained crowns; custom abutments require OEM components. Inclusive is a registered trademark of PrismaTik Dentalcraft, Inc. Glidewell HT and Hahn are trademarks of PrismaTik Dentalcraft, Inc. All other trademarks are property of their respective owners.