



4141 MacArthur Blvd. • Newport Beach, CA 92660
800-854-7256 • Fax 800-411-9722 • glidewell.com

Dr. Name _____ Phone # _____

Acct. # _____ Patient ID/Name _____

First Last

Address/Email _____ Deliver by 5 p.m. on _____

Enclosed with Case: Impressions Models Bite Photos Other: _____

PROVISIONAL RESTORATIONS
BioTemps Provisionals
Reinforcement: None Wire* Fiber
Transition C&B Smile Transitions
Abutment #(s) _____
Pontic #(s) _____ Total units _____
Splinted* Cement-On Implant
Individual Units Screw-Retained Implant
Amount of prep reduction: 1 mm* 2 mm
Perio treatment: Prepare tooth below gingival
on tooth #(s) _____ by _____ mm
Pontic site healing: Prepare ovate socket on tooth
#(s) _____ by _____ mm

DENTURES/FLIPPERS/FLEXIBLE PARTIALS
Denture Partial Select Phase
Handcrafted Flipper Custom tray
Digital (3D-printed) Valplast Bite rim
Immediate tcs Setup try-in
3D-printed DuraFlex Finish
Handcrafted
Copy (3D-printed)
Select Teeth
Digital teeth shade _____ Mould _____
Kenson Teeth
Standard for partials and handcrafted dentures
Shade _____ Mould _____
Premium Brand Teeth
extra charge for partials and handcrafted dentures
Shade _____ Brand _____ Mould _____

SIMPLY NATURAL METAL PARTIALS
Metal frame with acrylic and Kenson teeth
Frame Material Phase
SLM-printed cobalt chrome frame Metal frame try-in
Printed frame try-in
Esthetic Clasp Material
extra charge applies
Valplast/SLM-printed cobalt chrome frame
Frame w/occlus. rim
Frame w/setup try-in
Finish
Scan/Save File
extra charge applies
Lab select complete design

SNORING/SLEEP APNEA APPLIANCES
Buy 1 Buy 2 and save Scan/Save File
Silent Nite* Silent Nite with Glidewell Hinge
OASYS Hinge Appliance EMA
dreamTAP TAP 3 TL

ZIRCONIA/ALL-CERAMIC RESTORATIONS
BruxZir Full-Strength* (> 1,000 MPa) IPS e.max
NEW! BruxZir Esthetic (870 MPa) Lava Crowns & Bridges
Obsidian All-Ceramic Bilayered Clinical Zirconia

VENEERS
Obsidian veneer* NEW! BruxZir Esthetic veneer
IPS e.max veneer Layered IPS e.max veneer

PFM
Obsidian to Non-Precious* Obsidian to White High Noble
Obsidian to White Noble Porcelain fused to OcclusalGold YHN

SCREW-RETAINED RESTORATIONS
BruxZir Full-Strength* (w/ Ti-Base) Obsidian to White Noble
NEW! BruxZir Esthetic (w/ Ti-Base) Obsidian to White High Noble
IPS e.max (w/ Ti-Base) Bilayered Clinical Zirconia (w/ Ti-Base)

CUSTOM ABUTMENTS
Titanium* Zirconia w/ Ti-Base
Gold-Tone Titanium BIOMET 3i Encode
Gold Alloy Prepare existing abutment

COMPOSITE RESTORATIONS
Composite Fiber Reinforcement

FULL-CAST RESTORATIONS
Noble-Cast 45 YN (40% Au) White Noble
Noble-Cast 60 YHN (57.5% Au)* White High Noble (40% Au)
Non-Precious Post & Core

NIGHTGUARDS/RETAINERS/MIGRAINE PREVENTION
Upper Lower
Buy 1 Buy 2 and save Scan/Save File
Comfort3D (3D-printed, hard) Comfort H/S* (clear, hard with soft reline) Soft nightguard (buy 2 savings not available)
CLEARsplint (self-adjusting, hard) NTI-tss Plus NTI OmniSplint
Clear-Lock Retainers: Qty of 2 Qty of 4 Qty of 6
Clear ortho Hawley Essix Retainer (1 tooth)


PLAYSAFE MOUTHGUARDS
Jr Lt Lt Pro Med* Hvy Hvy Pro
Helmet strap _____ Specify color(s) on Rx
Name: _____

FINAL SHADE STUMP SHADE OCCLUSAL STAINING
Must indicate prepped tooth shade for all-ceramics
None Light*
Medium Dark

PONTIC DESIGN MARGIN AND METAL DESIGN
Labial Butt Junction Junct.Collar Excl.Cusp Incl.Cusp

IF NO OCCLUSAL CLEARANCE
Call doctor Spot opposing Metal occlusion Metal island Make this a permanent note

Rx Implant System _____ (If applicable)
Implant Diameter _____ mm
(See reverse for flat-rate pricing info)



Signature _____

License _____ Date _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. *Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.* Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewell.com/policies-and-warranties.



• BruxZir Restorations



• Custom Abutments



- All-Ceramic Restorations
- PFM Restorations
- Full-Cast Restorations




- Mandibular Advancement Devices
- Transition Crowns and Bridges®



- BioTemps® Provisionals
- Bite Splints
- Mouthguards and Nightguards

All rush cases must be prescheduled by calling 800-944-7874 before the case is shipped. Time of pickup and delivery may affect turnaround time.

PREPARATION GUIDELINES



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™
Certain®

CAMLOG®
SCREW-LINE

DENTSPLY Implants
ANKYLOS® C/X
ASTRA TECH Implant System®
ASTRA TECH Implant System® EV

Glidewell Direct
Hahn™ Tapered Implant System
Inclusive® Tapered Implant System

HIOSSEN®
HG System

MegaGen
AnyRidge® Implant System

Nobel Biocare
Brånemark System® RP
NobelActive®
NobelReplace®

Straumann®
Bone Level
Tissue Level

Zimmer Dental
Screw-Vent®

Restorations for the implant systems above are fabricated using components manufactured by Prisma™ Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. Inclusive is a registered trademark of Prisma™ Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prisma™ Dentalcraft, Inc. All other trademarks are property of their respective owners.