UNIVERSAL RX	

□ BioTemps Provisionals

☐ Transition C&B

on tooth #(s)

Abutment #(s)

Pontic #(s)

□ Splinted* ☐ Individual Units

Denture

☐ Handcrafted

☐ Immediate

☐ Digital (3D-printed)

O3D-printed

O Handcrafted

☐ Copy (3D-printed) Select Teeth Digital teeth shade

☐ Kenson Teeth

Shade

Frame Material

Esthetic Clasp Material

(extra charge applies)

□ Valplast/SLM-printed

□ tcs/SLM-printed

© 2024 Glidewell

cobalt chrome frame

cobalt chrome frame ☐ Lab select complete design

☐ Premium Brand Teeth



4141 MacArthur Blvd. • Newport Beach, CA 92660 800-854-7256 • Fax 800-411-9722 • glidewell.com

PROVISIONAL RESTORATIONS

Amount of prep reduction: □ 1 mm* □ 2 mm

☐ Pontic site healing: Prepare ovate socket on tooth

DENTURES/FLIPPERS/FLEXIBLE PARTIALS

Partial

☐ tcs

(Standard for partials and handcrafted dentures)

(extra charge for partials and handcrafted dentures) Brand

SIMPLY NATURAL METAL PARTIALS

Metal frame with acrylic and Kenson teeth

□ SLM-printed cobalt chrome frame □ Metal frame try-in

☐ Flipper

Valplast

□ DuraFlex

Mould

☐ Perio treatment: Prepare tooth below gingival

☐ Cast-Metal

□ Smile Transitions

☐ Cement-On Implant

☐ Screw-Retained Implant

by

Mould

Phase

☐ Finish

☐ Printed frame try-in

☐ Frame w/occlus. rim

☐ Frame w/setup try-in

(extra charge applies)

☐ Scan/Save File

Total units

Select Phase

☐ Bite rim

☐ Finish

☐ Custom tray

☐ Setup try-in

Dr. Name	Phone #
Acct. #	Patient ID/Name
Address/Email	····· -
Enclosed with Case: Impressions Mod	lels 🗆 Bite 🗅 Photos 🗅 Other:
ZIRCONIA/ALL-CERAMIC RESTORATIONS □ BruxZir Full-Strength* (> 1,000 MPa) □ NEW! BruxZir Esthetic (870 MPa) □ Obsidian All-Ceramic □ BruxZir Esthetic (870 MPa) □ Lava Crowns & Bridges □ Bilayered Clinical Zirconia	FINAL SHADE STUMP SHADE Must indicate prepped tooth shade for all-ceramics None Light* Medium Dark
VENEERS □ Obsidian veneer* □ IPS e.max veneer □ Layered IPS e.max veneer	PONTIC DESIGN MARGIN AND METAL DESIGN
PFM ☐ Obsidian to Non-Precious* ☐ Obsidian to White Noble ☐ Porcelain fused to OcclusalGold YHN	Labial Butt Junction Junct.Collar Excl.Cusp Incl.Cusp
SCREW-RETAINED RESTORATIONS BruxZir Full-Strength* Obsidian to White Noble (w/ Ti-Base) Obsidian to White High Noble NEW! BruxZir Esthetic Bilayered Clinical Zirconia	IF NO OCCLUSAL CLEARANCE □ Call doctor □ Spot opposing □ Metal occlusion □ Metal island □ Make this a permanent note
(w/ Ti-Base) (w/ Ti-Base) □ IPS e.max (w/ Ti-Base)	R Implant System
CUSTOM ABUTMENTS □ Titanium* □ Zirconia w/ Ti-Base □ Gold-Tone Titanium □ BIOMET 3i Encode □ Gold Alloy □ Prepare existing abutment	(If applicable) Implant Diameter mm (See reverse for flat-rate pricing info) mm 4 12 13 14
COMPOSITE RESTORATIONS □ Composite □ Fiber Reinforcement	2 15 ST 1 16 S
FULL-CAST RESTORATIONS □ Noble-Cast 45 YN (40% Au) □ White Noble □ Noble-Cast 60 YHN (57.5% Au)* □ White High Noble (40% Au) □ Non-Precious □ Post & Core	32 17 31 18 53 30 26 23 19 54 1
NIGHTGUARDS/RETAINERS/MIGRAINE PREVENTION Upper Lower Buy 1 Buy 2 and save Scan/Save File	29 25 24 20 21 27 22 27 27
□ Comfort3D □ Comfort H/S* □ Soft nightguard (clear, hard with soft reline) □ Comfort H/S*	
□ CLEARsplint □ NTI-tss Plus □ NTI OmniSplint (self-adjusting, hard)	
☐ Clear-Lock Retainers: ○ Qty of 2 ○ Qty of 4 ○ Qty of 6	
☐ Clear ortho ☐ Hawley ☐ Essix Retainer (1 tooth)	Signature
PLAYSAFE MOUTHGUARDS	License Date
□ Jr □ Lt □ Lt Pro □ Med* □ Hvy □ Hvy Pro	
☐ Helmet strapSpecify color(s) on Rx Name:	Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details

SNORING/SLEEP APNEA APPLIANCES

O Buy 1	O Buy 2 a	nd save	O Scan/Save F	-ile
☐ Silent Nite*		☐ Silent	Nite with Glidewel	l Hing
OASYS Hinge	Appliance	□ EMA	☐ flexTAP	

☐ TAP 3 TL □ dreamTAP

*Standard unless specified otherwise

Name:

details.

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewell.com/policies-and-warranties.







- All-Ceramic Restorations
- PFM Restorations
- Full-Cast Restorations



- BioTemps® **Provisionals**
- Bite Splints
- Mouthguards and Nightguards



Silent Nite

flexTAP

- dreamTAP
- TAP 3 TL

Silent Nite with


Glidewell Hinge

EMA

All rush cases must be prescheduled by calling 800-944-7874 before the case is shipped. Time of pickup and delivery may affect turnaround time.

PREPARATION GUIDELINES

Anterior LINGUAL LABIAL



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™ Certain®

CAMLOG® SCREW-LINE **DENTSPLY Implants** ANKYLOS® C/X

ASTRA TECH Implant System® ASTRA TECH Implant System® EV

> **Nobel Biocare** Brånemark System® RP NobelActive® NobelReplace®

Glidewell Direct

Glidewell HT™ Implant System Hahn[™] Tapered Implant System Inclusive® Tapered Implant System

Straumann® **Bone Level** Tissue Level

HIOSSEN® HG System

Zimmer Dental Screw-Vent®

MegaGen AnyRidge® Implant System

Restorations for the implant systems above are fabricated using components manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems above are fabricated using components manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems above are fabricated using components manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems above are fabricated using components manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems above are fabricated using components manufactured by Prismatik Dentalcraft, Inc. Prices may vary for other implant systems are fabricated using components. components. Inclusive is a registered trademark of Prismatik Dentalcraft, Inc. Glidewell HT and Hahn are trademarks of Prismatik Dentalcraft, Inc. All other trademarks are property of their respective owners.