

WHEN REALISM MATTERS MOST



Prescribe for
\$149*
per unit

*Jim Glidewell, CDT
Founder and President, Glidewell
BruxZir Glass Patient*

"BruxZir® Glass isn't just an advancement, it's the culmination of everything I've worked for in my half-century career — a crown that's virtually indistinguishable from natural teeth. Lifelike beauty, and all the strength you'll ever need."

- **A Better Way to Make a Crown** — Enamel-on-dentin design, with lithium disilicate pressed over zirconia
- **Show Only What You Want** — Wraparound translucency replicates the vitality of real teeth
- **And Nothing You Don't** — A strong, anatomically shaped core masks even the darkest stumps


*Price does not include shipping or applicable taxes.

All BruxZir Glass photos are unedited and unenhanced.
Crowns #6–11 were crafted by Thang Tran, CDT.




**Prescribe
BruxZir® Glass**
\$149 Per Unit

 **Glidewell**
for the sake of smiles



Crowns That Look Like Teeth, Not Crowns

BruxZir® Glass is built for realism. With nature as inspiration, this new restorative solution uses an enamel-on-dentin design — with a lithium disilicate “enamel” layer pressed over a BruxZir Zirconia “dentin” layer.



Esthetic Advantages of the Glass-on-BruxZir Design

- Anatomical depth — with internal morphology that bleeds through the outer surface
- True-to-life translucency in the incisal and interproximal
- Mask even the darkest stumps and discolorations, without sacrificing the final look



**Prescribe
BruxZir® Glass**
\$149 Per Unit

 **Glidewell**
for the sake of smiles