



IMPLANT SOLUTIONS

"My family has used Glidewell for over 5 decades. Glidewell has always been at the forefront of implant dentistry, making it easier for doctors to take their implant practices to a very high level of success and predictability."

—Randolph R. Resnik, DMD, MDS

FULL-ARCH IMPLANT PROSTHESIS



\$2,795 / \$1,398

SCREW-RETAINED CROWN



\$339 / \$170

**50%
OFF**

CUSTOM ABUTMENT



\$318 / \$159

GlidePRO™



\$2,750 / \$1,375

SPECIAL PRICING: SAVE 50% ON YOUR FIRST CASE

glidewell.com | 800-839-9755

SPECIAL SAVINGS FOR NEW CUSTOMERS ON OUR MOST POPULAR PRODUCTS



GlidePRO™

- All-on-X stackable surgical guide system
- Same-day design of immediate provisional for “teeth-in-a-day”
- Print in-office or have Glidewell overnight a printed provisional

~~\$2,750/arch list price~~ | **\$1,375** Tier 2



DIGITAL BRUXZIR® FULL-ARCH IMPLANT PROSTHESIS

- 100% monolithic zirconia, backed by lifetime warranty
- Stability and function of a fixed implant restoration
- Includes provisional try-in implant prosthesis

~~\$2,795/arch list price~~ | **\$1,398/arch**[†]



DIGITAL BRUXZIR ESTHETIC FULL-ARCH IMPLANT PROSTHESIS

- Enhanced translucency for natural-looking, lifelike esthetics
- 100% monolithic zirconia, backed by lifetime warranty
- Includes provisional try-in implant prosthesis

~~\$2,995/arch list price~~ | **\$1,498/arch**[†]

Promo available only to first-time customers to product.

*Price does not include round-trip shipping or applicable taxes. See reverse of Rx for specific implant systems eligible for flat-rate pricing. †Flat-rate lab fee available for most major implant systems.



4141 MacArthur Blvd. • Newport Beach, CA 92660
 866-861-7955 • Fax 800-411-9722 • glidewell.com

Dr. Name _____
 Phone # _____ Acct. # _____
 Patient ID/Name _____
First Last
 Address/Email _____

Deliver by 5 p.m. on _____ Call before starting case

Enclosed with case: Impressions Models Bite Photos Other: _____

Rx Implant Information

Tooth #	Diameter	Implant System
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____

\$80 off Mini Implant Overdenture



\$100 off Locator Overdenture



Upper Lower **Tooth Shade** _____
Gingival Shade _____

Stage of Service Needed:

- Wax rim Implant verification jig (IVJ)
- Wax setup Reset
- Provisional implant prosthesis
- Final BruxZir Esthetic Full-Arch Implant Prosthesis (see reverse for flat-rate pricing info and limited warranty details)
- Final BruxZir Full-Arch Implant Prosthesis (see reverse for flat-rate pricing info and limited warranty details)

20% off Screw-Retained Hybrid Denture



50% off BruxZir Full-Arch Implant Prosthesis



50% off GlidePRO Surgical Guides



50% off BruxZir Esthetic Full-Arch Implant Prosthesis



Signature _____

License # _____ Date _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

[†]Flat-rate pricing is available for most major implant systems (see reverse) and covers all appointments, parts, materials and labor.

SPECIAL OFFER*

- Carefully package your case, including this Rx, and tape box securely closed **or submit your case digitally with the new Digital Case Submission on My Account.**
- To schedule shipping pickup, call us at **800-854-7256.**
- Discounts can be applied toward two arches per Rx.

Offer expires December 31, 2026.

**Price does not include round-trip shipping or applicable taxes. No copies accepted. Offer applicable for a maximum of one specially priced Rx per account for specified products only. Special pricing not valid with any other offer and cannot be combined with reduced-pricing offer for restorations over Glidewell HT Implants.*

BRUXZIR IMPLANT PROSTHESIS

- BruxZir Esthetic Full-Arch Implant Prosthesis[†]**
(Lifetime warranty)
 (Precision-milled esthetic solid zirconia with tooth and gingival tissue shade. NOTE Complete service includes provisional try-in implant prosthesis; a duplicate provisional can be purchased for an extra fee.)
- BruxZir Full-Arch Implant Prosthesis (Full-Strength)[†]**
(Lifetime warranty)
 (Precision-milled solid zirconia with tooth and gingival tissue shade. NOTE: Complete service includes provisional try-in implant prosthesis; a duplicate provisional can be purchased for an extra fee.)
- Provisional Try-In Implant Prosthesis**
(Provisional with tooth and gingival tissue shade)

Tooth Setup:

- Ideal Characterized
- Copy study model Copy existing denture
- Add lip support
- Male Female Age _____

Gingival Shade:

- G00 (Lightest) G0 (Light)
- G1 (Standard) G3 (Medium)
- G4 (Dark) G5 (Darkest)

Bite splint (Additional fee applies)

IMPLANT DENTURE APPLIANCE

- Locator Overdenture[†]**
(Kenson teeth standard)
- Mini Implant Overdenture**
(Kenson teeth standard)
- Screw-Retained Hybrid Denture[†]**
(Premium teeth standard)
- Kenson teeth** (Included at no extra charge)
Shade _____
Mould _____
- Premium brand teeth** (Extra charge may apply)
Shade _____ **Brand** _____
Mould _____
- Name on appliance (Additional charge)

DENTURE REINFORCEMENT

(overdenture cases only)

- Horseshoe cast palate**
- Cast MESH 360° wraparound**

TISSUE SURFACE DESIGN



- Convex[†]
- Cover exposed implant



lingual

- Modified convex
- Provide floss space

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewelldental.com/policies-and-warranties.



• Glidewell HT™
Implants



• Custom
Abutments



• BruxZir®
Restorations

All rush cases must be prescheduled by calling **800-944-7874** before the case is shipped. Time of pickup and delivery may affect turnaround time.

IMPLANT WORKING TIMES

Please allow full working time for each product selected. If case includes both upper and lower arches, please allow double working time for case completion. Working times are not guaranteed and do NOT include weekends or holidays.

BruxZir® Full-Arch Implant Prosthesis

Not all steps are required for Digital Workflows.

Days in Lab

Implant verification jig, custom tray and wax rim	9
Wax setup.....	5
Provisional try-in implant prosthesis	6
Final prosthesis	8

FLAT-RATE PRICE ON THE BRUXZIR® IMPLANT PROSTHESIS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

**BIOMET 3i™
Certain®**

**CAMLOG®
SCREW-LINE**

**DENTSPLY Implants
ANKYLOS® C/X
ASTRA TECH Implant System®
ASTRA TECH Implant System® EV**

**Glidewell Direct
Glidewell HT™ Implant System
Hahn™ Tapered Implant System
Inclusive® Tapered Implant System**

**HIOSSEN®
HG System**

**MegaGen
AnyRidge® Implant System**

**MIS®
C1 Implant System**

**Nobel Biocare
Brånemark System® RP
NobelActive®
NobelReplace®**

**Straumann®
Bone Level
Bone Level SC
Tissue Level**

**Zimmer Dental
Screw-Vent®**

Restorations for the implant systems above are fabricated using components manufactured by Prismatic Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. Inclusive is a registered trademark of Prismatic Dentalcraft, Inc. Glidewell HT and Hahn are trademarks of Prismatic Dentalcraft, Inc. All other trademarks are property of their respective owners.



4141 MacArthur Blvd. • Newport Beach, CA 92660
 866-861-7955 • Fax 800-411-9722 • glidewell.com

SPECIAL OFFER*

- Carefully package your case, including this Rx, and tape box securely closed.
- To schedule shipping pickup, call us at **800-854-7256** or log in to myaccount.glidewell.com.
- Use this Rx to get 50% off screw-retained crowns, implant crowns and custom abutments, or screwmentable restorations (limit 4 units per Rx).

Offer expires December 31, 2026. Price does not include round-trip shipping or applicable taxes.
 Limit one special offer per implant case. No copies accepted. Special pricing not valid with any other offer.

Dr. Name _____ Acct. # _____
 Phone # _____ Email _____
 Address _____ City/State/ZIP _____
 Patient ID/Name _____ First _____ Last _____
Deliver by 5 p.m. on _____ Call before starting case

SHADE INSTRUCTIONS

Tooth # _____ Final Shade _____
 Occlusal Staining: Light** Med
 Dark None

RESTORATION PONTIC DESIGN




SELECT CUSTOM ABUTMENT

Titanium** Zirconia w/Ti-Base
 Gold-Colored Titanium

PARALLEL ABUTMENTS

No Yes (indicate which abutments will have restorations splinted together for insertion)

ABUTMENT MARGIN DEPTH



If left blank, default values will be used

ABUTMENT EMERGENCE PROFILE



Surgical Tissue No tissue
 placement displacement** displacement

**Standard unless specified otherwise

IF NO OCCLUSAL CLEARANCE

Call doctor Spot opposing
 Make this a permanent note in my master file

Rx



Implant Information

Tooth #	Diameter	Implant System
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____
_____	_____ mm	_____



(See reverse for flat-rate pricing info)

SELECT RESTORATION TYPE

Crown Splinted Crowns Bridge

CEMENT-RETAINED RESTORATION

BruxZir Full-Strength** (> 1,000 MPa)
 BruxZir Esthetic (870 MPa)
 NEW! BruxZir Fusion
 IPS e.max

SCREW-RETAINED RESTORATION

BruxZir Full-Strength** (w/ Ti-Base)
 BruxZir Esthetic (w/ Ti-Base)
 NEW! BruxZir Fusion (w/ Ti-Base)
 IPS e.max (w/ Ti-Base)

SCREWMENTABLE

(cementable crown with screw-access hole over custom abutment)

BruxZir Full-Strength**
 BruxZir Esthetic
 NEW! BruxZir Fusion
 IPS e.max

CONTOUR AND OCCLUSION DESIGN

Embrasures: Closed** Open
 Occlusion: Light** Ideal
 Open _____ mm Out
 Contacts: Broad & Tight** Pinpoint
 Light



50% off screw-retained crowns, implant crowns and custom abutments, or screwmentable restorations

Signature _____

License # _____ Date _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewell.com/policies-and-warranties.



- Glidewell HT™ Implants



- Custom Abutments



- BruxZir® Restorations



- All-Ceramic Restorations
- PFM Restorations

All rush cases must be prescheduled by calling 800-944-7874 before the case is shipped. Time of pickup and delivery may affect turnaround time.

FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

**BioHorizons®
Tapered Internal**

**BIOMET 3i™
Certain®**

**CAMLOG®
CONELOG®
SCREW-LINE**

**Dentium®
NR Line**

**DENTSPLY Implants
ANKYLOS® C/X
ASTRA TECH Implant System®
ASTRA TECH Implant System® EV**

**Glidewell Direct
Glidewell HT™ Implant System
Hahn™ Tapered Implant System
Inclusive® Tapered Implant System**

**HIOSSEN®
HG System**

**Implant Direct®
InterActive®
Legacy®
RePlant®**

**Keystone® Dental
PrimaConnex®**

**MegaGen
AnyRidge® Implant System**

**MIS®
C1 Implant System
SEVEN®**

**NeoDent®
Grand Morse® Implant System**

**Neoss®†
ProActive® Narrow
ProActive Standard**

**Nobel Biocare
Brånemark System® RP
NobelActive®
NobelReplace®**

**Straumann®
BLX
Bone Level
Bone Level SC
Tissue Level**

**Zimmer Dental
Eztetic®
Screw-Vent®**

Restorations for the implant systems above are fabricated using components manufactured by PrismaTik Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. †For Neoss®, flat-rate pricing is only available for screw-retained crowns; custom abutments require OEM components. Inclusive is a registered trademark of PrismaTik Dentalcraft, Inc. Glidewell HT and Hahn are trademarks of PrismaTik Dentalcraft, Inc.

All other trademarks are property of their respective owners.

THE STRONGEST WARRANTY IN IMPLANT DENTISTRY



CUSTOM ABUTMENTS

- Optimizes emergence profile and soft-tissue contours
- Available in titanium and zirconia with titanium base
- Eases the removal of excess cement

~~\$318/unit list price~~ | **\$159/unit*** (up to 4 units)



SCREW-RETAINED RESTORATIONS

- Eliminate cement problems around implants
- Facilitate natural emergence profile and gingival margins
- Seat quickly and precisely so you can move on to your next patient

~~\$339/unit list price~~ | **\$170/unit*** (up to 4 units)



SCREWMENTABLE RESTORATION

- Benefits of custom abutment without cement clean up; cemented in lab
- Retrievability of the implant restoration
- Lifetime warranty of abutment and crown

~~\$475/unit list price~~ | **\$238/unit*** (up to 4 units)

Promo available only to first-time customers to product.

INTRODUCING

Glidewell HT™
IMPLANT SYSTEM

THE GLIDEWELL HT™

NO-CHARGE

IMPLANT PROGRAM



Restore a Glidewell HT Case,
Earn a No-Charge Implant Credit

SCAN TO LEARN MORE



Glidewell HT™
IMPLANT SYSTEM

Official implant of the



*Price does not include round-trip shipping or applicable taxes. *Flat-rate lab fee available for most major implant systems. See reverse of Rx for specific implant systems eligible for flat-rate pricing.



GL-6700081-030126