

# The combination of speed, esthetics and durability you desire

- Recommended for inlays, onlays, and anterior and posterior crowns
- Mills in under 10 minutes\*
- Available exclusively for Glidewell milling systems
- Long-lasting shade and gloss retention



RETAINS SHADE



RESISTS STAINS



NO OVEN NEEDED

Before



This patient presented with failing direct composite restorations on teeth #30 and #31. A CAMouflage® NOW A1 overlay on tooth #30 and MOBL onlay on tooth #31, provided a minimally invasive single-visit solution.

- CAMouflage NOW is an excellent restorative alternative for indirect dentistry. I can provide high-quality conservative restorations to my patients in a matter of minutes. This, in turn, has increased the productivity of my practice.
- Chad Duplantis, DDS | Fort Worth, TX



# **CAMouflage® NOW**

#### **BENEFITS**

- Shade stability long after delivery
- Rapid patient delivery with no oven wait
- Choose from a wide range of lifelike shades that blend well with existing dentition



Sold in 5-block packs

# **READY TO DELIVER IN UNDER 10 MINUTES**



**MILL** 



NO FINISHING/ POLISHING REQUIRED



NO OVEN TIME REQUIRED



**TOTAL TIME** 

### **CIFICATIONS**

#### **MATERIAL**

- Approximately 80% ceramic-filled nanohybrid resin composite
- Flexural strength: 193 MPa

#### **IDEAL FOR**


- Single-unit anterior and posterior crowns
- Inlays and onlays

#### **SHADES**

- 7 VITA® Classical Shades, 5 in 2 translucencies: A1HT/LT, A2HT/LT, A3HT/LT, A3.5HT/LT, B1HT/LT, C2LT, D2LT
- 1 VITA Bleach Shade: BLLT

# **EASY PREPARATION**

#### **ONLAY**





- Prepare cavity walls tapered 5-6° to the long axis.
- Wall thickness should be ≥1.5 mm.
- Isthmus thickness should be ≥1.5 mm.

\*Milling times vary based on prescription and mill used. Average milling time is shown. VITA is a registered trademark of VITA Zahnfabrik.

MKT-10487\_6 PK-280353-032322

