

BruxZir[®] NOW SRC



MAXIMUM
STRENGTH



NO OVEN
NEEDED



PRE-CEMENTED
ABUTMENT



< 60 MIN
MILL TIME*

THE FULLY SINTERED SCREW-RETAINED ZIRCONIA IMPLANT CROWN FOR IN-OFFICE USE

BENEFITS

- Deliver implant crowns right from the mill in a single appointment
- Pre-cemented titanium abutment allows for easy delivery
- CAD/CAM fabrication ensures esthetic gingival margins and a natural emergence profile

FOR MORE INFORMATION

glidewell.com/io-materials | 844-949-7184

glidewell.io[™]
IN-OFFICE SOLUTION



This BruxZir® NOW SRC was digitally designed and milled chairside, achieving ideal contacts, occlusion, and emergence profile without the need for additional laboratory collaboration.

SPECIFICATIONS

MATERIAL

- Pre-shaded, fully sintered zirconia
- Pre-installed titanium base
- Flexural strength: >1,000 MPa

IDEAL FOR


- Prosthetic rehabilitation

SHADES

- 7 VITA Classical Shades (A1, A2, A3, A3.5, B1, C2, D2)

INSTRUCTIONS FOR USE

- 1 Fully sintered zirconia superstructure
- 2 Central screw access hole
- 3 Pre-cemented titanium abutment
- 4 Plastic mandrel
- 5 Abutment screw



IMPLANT SYSTEM	IMPLANT Ø	PLATFORM Ø
Glidewell HT™ Implant System	Ø3.5 mm	Ø3.5/4.3 mm
	Ø4.3 mm	Ø3.5/4.3 mm
	Ø5.0 mm	Ø5.0 mm
	Ø7.0 mm	Ø7.0 mm
Inclusive Tapered Implant System	Ø3.2 mm	Ø3.0 mm
	Ø3.7 mm	Ø3.5 mm
	Ø4.2 mm	Ø3.5 mm
	Ø4.7 mm	Ø4.5 mm
	Ø5.2 mm	Ø4.5 mm

FOR MORE INFORMATION

glidewell.com/io-materials | 844-949-7184

*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill. VITA is a registered trademark of VITA Zahnfabrik. The Glidewell HT Implant System is manufactured by PrismaDentalcraft, Inc., a wholly owned subsidiary of Glidewell Laboratories. Glidewell HT is a trademark of PrismaDentalcraft, Inc.

PK-5158208-061926