



With the glidewell.io™ In-Office Solution, effortlessly mill the crowns you want, when you want, with these special features:

- Mill fully sintered BruxZir® crowns and bridges with no oven time required — just polish and cement, and you're done.
- Let MarginAI™ and CrownAI™ do the design work for you so you can spend your time with your patients.
- You're connected directly to Glidewell for training and support at any time.
- Love your scanner? Great! — glidewell.io is an open system and works with your preferred scanning system.

“The glidewell.io In-Office Solution brings the ‘wow’ factor. In just one visit, your patient is receiving a beautiful BruxZir NOW crown that’s going to last a lifetime.”

– Christina Do, DDS
Costa Mesa, CA



“Who better than the world’s largest dental lab that makes these materials day in and day out to know what the materials should look like at the end of the day? The service was there. The material was there. The organization was backed. For us, it was very logical and just made a ton of sense.”

– Oussama Founas, DDS
Chicago, IL



“The glidewell.io system and its incorporated software are AMAZING! Scanning the patient and designing the restoration are easy to learn. The blocks mill quickly and come out of the mill basically polished and ready to go — CEREC® could never match that!”

– Robert D. Rutner, DDS, DICOI
Mountain View, CA



FASTMILL.IO™ PACKAGE

Includes:

- fastmill.io™ In-Office Mill
- fastdesign.io™ Software and Design Station
- Mill Accessory Kit
- One- or five-year Standard Warranty
- One- or five-year fastdesign.io Software License
- 45 Free Milling Blocks
- Free glidewell.io™ Level I Advanced Course at CA, NC, TX or FL

MILLING MATERIALS

The glidewell.io™ In-Office Solution offers dentists a versatile range of restorative materials designed to meet the specific needs of any clinical situation, whether for a full-coverage crown or a bridge. The fastmill.io™ In-Office Mill stands out as the only machine capable of producing same-day BruxZir® NOW—our fully sintered, oven-free restoration material.



ZIRCONIA MATERIAL

High strength • Posterior confidence • Implant-ready options

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
BruxZir® NOW		✓						> 1,000 MPa	< 35 minutes (single-unit crowns)
BruxZir® NOW BRIDGE BLOCK			✓					> 1,000 MPa	< 2.5 hours (single-unit crowns)
BruxZir® esthetic NOW	✓	✓		✓				870 MPa	< 35 minutes
BruxZir® esthetic NOW POSTERIOR		✓						870 MPa	< 40 minutes
BruxZir® NOW SRC		✓				✓		> 1,000 MPa	< 40 minutes

RESIN & PMMA

Fast • Versatile • Ideal for temporaries and conservative cases

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
CAMouflage® NOW	✓	✓		✓	✓			193 MPa	< 10 minutes
BioTemps® NOW	✓	✓		✓	✓		✓	≥ 50 MPa	< 10 minutes
BioTemps® NOW BRIDGE BLOCK			✓				✓	≥ 50 MPa	< 30 minutes (3-unit bridges)

GLASS CERAMIC MATERIALS

Esthetics-driven • Conservative • Natural translucency

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
Obsidian® MILLING BLOCKS	✓	✓		✓	✓			385 MPa	< 20 minutes
Obsidian® NOW	✓							476 MPa	< 35 minutes



LEARN MORE

glidewell.com/fastmillio-in-office-mill
844-949-7184

LEARN MORE ABOUT OUR BRUXZIR CEMENT

glidewelldirect.com/products/bruxzir-dual-cure-resin-cement

