

BioTemps® **NOW**

SCREW-RETAINED



RETAINS
GINGIVAL POSITION



NO OVEN
NEEDED



PATIENT-PLEASING
RESULTS



< 15 MIN
MILL TIME*

PMMA MATERIAL FOR IN-OFFICE SCREW-RETAINED PROVISIONALS

BENEFITS

- Provide patients with a same-day preview of their final restoration
- Support soft tissue contouring and gingival tissue management during healing
- Reduce chairtime with a retrievable screw-retained provisional for easy adjustment and evaluation

FOR MORE INFORMATION

glidewell.com/io-materials | 844-949-7184

glidewell.io™
IN-OFFICE SOLUTION



This patient presented with failing restorations and recurrent decay on teeth #12–15. Because of the number of units involved, BioTemps® NOW provisional restorations were milled and seated for a durable and comfortable solution while awaiting final treatment.

SPECIFICATIONS

MATERIAL

- PMMA
- Flexural strength:
≥ 50 MPa

IDEAL FOR

- Temporary full-coverage anterior and posterior crowns

SHADES

- 6 VITA® Classical Shades, 1 Bleach Shade: A1, A2, A3, A3.5, B1, C2, BL

IMPLANT SYSTEM	SIZES	IMPLANT Ø	PLATFORM Ø
Glidewell HT™ Implant System		Ø3.5 mm Ø4.3 mm Ø5.0 mm Ø7.0 mm	Ø3.5/4.3 mm Ø3.5/4.3 mm Ø5.0 mm Ø7.0 mm
Inclusive Tapered Implant System		Ø3.2 mm Ø3.7 mm Ø4.2 mm Ø4.7 mm Ø5.2 mm	Ø3.0 mm Ø3.5 mm Ø3.5 mm Ø4.5 mm Ø4.5 mm
Nobel Biocare NobelActive®	NP	Ø3.5 mm	Ø3.5mm
	RP	Ø4.3 mm	Ø4.3mm
		Ø5.0 mm	
Zimmer Dental Screw-Vent®	3.5	Ø3.7 mm	Ø3.5 mm
		Ø4.1 mm	
	4.5	Ø4.7 mm	Ø4.5 mm
Straumann® Bone Level	NC	Ø3.3 mm	Ø3.3mm
	RC	Ø4.1 mm	Ø4.1mm
		Ø4.8 mm	
Hiossen® HG System	Mini	Ø3.5 mm	Ø3.5 mm
	Standard/Regular	Ø4.0 mm	Ø4.0mm
		Ø4.5 mm	
		Ø5.0 mm	
		Ø5.5 mm	

FOR MORE INFORMATION

glidewell.com/io-materials | 844-949-7184

*Milling times vary based on prescription and mill used. Average milling time is shown. VITA is a registered trademark of VITA Zahnfabrik. NobelActive is a registered trademark of the Nobel Biocare group. Screw-Vent is a registered trademark of Zimmer Dental Inc. Straumann is a registered trademark of Straumann Holding AG. HIOSEN is a registered trademark of OSSTEM IMPLANT CO., LTD.

PK-7658000-100126