

SINGLE-VISIT DENTISTRY WITH ZERO HASSLE

Learn how the glidewell.io™ In-Office Solution can improve your workflow, boost patient satisfaction and deliver an impressive ROI for your practice.



LEARN MORE

glidewell.io/fastmillio-in-office-mill
844-949-7184

glidewell.io™
IN-OFFICE SOLUTION

The glidewell.io™ In-Office Solution is a complete chairside CAD/CAM system that makes it easy to design, mill and deliver same-day crowns.

- Mill fully sintered BruxZir® crowns and bridges with no oven time required — just polish and cement, and you're done.
- Let MarginAI™ and CrownAI™ design restorations for you so you don't have to do the role of a technician.
- Get free training and support from a Glidewell technician at any time.
- Love your scanner? Great! — glidewell.io is an open system and works with your preferred intraoral scanner.



"The ROI has been significant. Because our lab bills are a fraction of what they used to be, fastmillio™ allows us to reinvest in our practice over and over again."

– Thomas Fuchs, DMD
Louisville, KY

"The support with Glidewell™ is amazing. Especially if you're just starting with digital, they will happily hold your hand and guide you through everything."

– Timothy Loo, DDS
San Francisco, CA



"What loved about glidewell.io™ was how it reduced the steps. You scan, design and mill - all in one streamlined system. The blocks are fully sintered zirconia, so there's no baking, glazing, or firing."

– Ousamma Founas, DDS
Chicago, IL

FASTMILL.IO™ PACKAGE

Includes:

- fastmill.io™ In-Office Mill
- fastdesign.io™ Software and Design Station
- Mill accessory kit
- One- or five-year standard warranty
- One- or five-year fastdesign.io software license
- 45 free milling blocks
- Free glidewell.io Level I Advanced Course at CA, NC, TX or FL

Opinions expressed in testimonials are solely the opinion of the doctor referenced and not those of Glidewell or its representatives.
CEREC is a registered trademark of Dentsply Sirona, Inc.

IDEAL RESTORATIVE MATERIALS FOR EVERY INDICATION

The glidewell.io™ In-Office Solution offers a range of restorative materials suited for virtually any clinical situation, from full-coverage crowns to bridges. As part of this integrated ecosystem, the fastmill.io™ In-Office Mill is the only machine capable of producing same-day BruxZir® NOW—our fully sintered, oven-free zirconia restoration material.

ZIRCONIA MATERIAL

High strength • Posterior confidence • Implant-ready options

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
BruxZir® NOW		✓						> 1,000 MPa	< 35 minutes (single-unit crowns)
BruxZir® NOW BRIDGE BLOCK			✓ (3-4 units)					> 1,000 MPa	< 2.5 hours (3-unit bridges)
BruxZir® esthetic NOW	✓	✓ (premolars only)		✓				870 MPa	< 35 minutes
BruxZir® esthetic NOW POSTERIOR		✓						870 MPa	< 40 minutes
BruxZir® NOW SRC		✓				✓		> 1,000 MPa	< 40 minutes

RESIN & PMMA

Fast • Versatile • Ideal for temporaries and conservative cases

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
CAMouflage® NOW	✓	✓		✓	✓			193 MPa	< 10 minutes
BioTemps® NOW	✓	✓		✓	✓		✓	≥ 50 MPa	< 10 minutes
BioTemps® NOW BRIDGE BLOCK			✓ (3-4 units)				✓	≥ 50 MPa	< 30 minutes (3-unit bridges)

GLASS CERAMIC MATERIALS

Esthetics-driven • Conservative • Natural translucency

Materials	Anterior Single Unit	Posterior Single Unit	Bridges	Veneers	Inlay/Onlay	Implant	Temporary	Flexural Strength	Milling Time
Obsidian® MILLING BLOCKS	✓	✓		✓	✓			385 MPa	< 20 minutes
Obsidian® NOW	✓							476 MPa	< 35 minutes



LEARN MORE ABOUT OUR MATERIALS

glidewell.com/io-materials

BruxZir® NOW



BruxZir® NOW BRIDGE BLOCK



BruxZir® esthetic NOW



BruxZir® esthetic NOW POSTERIOR



CAMouflage® NOW



BioTemps® NOW



BioTemps® NOW BRIDGE BLOCK




Obsidian® MILLING BLOCKS



Obsidian® NOW



BruxZir® NOW SRC



Clinical dentistry by Justin Chi, DDS, CDT and Jinny Bender, DMD



LEARN MORE

glidewell.io/fastmillio-in-office-mill
844-949-7184

LEARN MORE ABOUT OUR BRUXZIR CEMENT

glidewelldirect.com/products/bruxzir-dual-cure-resin-cement

