

**BRUXZIR RX**



22027 70th Ave. S • Kent, WA 98032  
833-814-7878 • Fax 978-313-6374 • [glidewell.com/local](http://glidewell.com/local)

**SPECIAL BRUXZIR ZIRCONIA OFFER**

1. Carefully package your case, including this Rx, and tape securely closed.
2. To schedule a pickup, call us at **833-814-7878**.
3. Please allow up to five working days in lab.

Dr. Name \_\_\_\_\_ Acct. # \_\_\_\_\_

Address \_\_\_\_\_ City/State/ZIP \_\_\_\_\_

Patient ID/Name \_\_\_\_\_ Age \_\_\_\_\_

First \_\_\_\_\_ Last \_\_\_\_\_


Deliver by 5 p.m. on \_\_\_\_\_  Call before starting case

Enclosed with case:  Impressions  Bite  Models  Articulator  Shade Tab  D-Wax  Pre-Op Models  Photos

**WEB Rx**



**NOTE:** Please send a study model on all work involving anterior teeth.




Signature \_\_\_\_\_ License # \_\_\_\_\_ Date \_\_\_\_\_

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

**CHOOSE PRODUCT**

- BruxZir Full-Strength\*\* (> 1,000 MPa)
- NEW!** BruxZir Esthetic (870 MPa)

**FINAL CERAMIC SHADE**



Indicate Shade Here \_\_\_\_\_

**PRESENT TOOTH OR STUMP SHADE**




Indicate Shade Here \_\_\_\_\_

**OCCUSAL STAINING**

- Light\*\*  Med  Dark  None

**INCISAL SHAPE INSTRUCTIONS**



- Rounded  Squared  Pointed

**PONTIC DESIGN**



\*\*Standard unless specified otherwise

## TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.




**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [glidewell.com/policies-and-warranties](http://glidewell.com/policies-and-warranties).

## PREPARATION GUIDELINES



### BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

### BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

## THE JIM GLIDEWELL SIGNATURE ANTERIOR SERIES

Choose the desired esthetic outcome for your patient

### SQUARE-TAPERED



Exemplify boldness and youthfulness.

### SOFT-SQUARED



Show energetic professionalism.

### OVOID



Convey charm and softness.

### TRAPEZOID



Demonstrate confidence.

### SQUARED



Strike an athletic tone.

### TRIANGLE-TAPERED



Exhibit mature experience.