



• IMPLANT FIXED Rx •

GLIDEWELL LABORATORIES

4141 MacArthur Blvd. • Newport Beach, CA 92660

800-839-9755 • Fax 800-411-9722
glidewell dental.com

Dr. Name _____ Phone # _____

Acct. # _____ Patient Name _____
First Last

Address/Email _____ Deliver by 5 p.m. on _____ See Reverse for Working Times

Enclosed with Case: Impressions Models Bite Photos Other: _____

Rx Indicate implant system _____

Indicate implant size _____ mm

Please order all implant components for this case



SHADE INSTRUCTIONS



Tooth No. _____

Final Shade _____

Occlusal Staining: Light* Med Dark None

SELECT FINAL CUSTOM ABUTMENT

- | | |
|--|---|
| Choose Material | Choose Brand |
| <input type="checkbox"/> Titanium* | <input type="checkbox"/> Inclusive* |
| <input type="checkbox"/> Gold-Colored Titanium | <input type="checkbox"/> BIOMET 3i Encode |
| <input type="checkbox"/> Zirconia w/ Ti-Base | |
| <input type="checkbox"/> Gold Alloy | |


Prepare existing abutment

See reverse for supported implant systems

PARALLEL ABUTMENTS


No Yes (indicate which abutments will have restorations splinted together for insertion)

ABUTMENT MARGIN DEPTH




If left blank, default values will be used

ABUTMENT MARGIN DESIGN



ABUTMENT EMERGENCE PROFILE



SELECT RESTORATION TYPE

Crown Splinted Crowns Bridge

CEMENT-RETAINED RESTORATIONS

- | | |
|---|---|
| <input type="checkbox"/> BruXZir Full-Strength* (1,150 MPa) | <input type="checkbox"/> IPS e.max (anterior) |
| <input type="checkbox"/> NEW! BruXZir Esthetic (870 MPa) | <input type="checkbox"/> IPS e.max (posterior) |
| <input type="checkbox"/> PrismaTik Clinical Zirconia | <input type="checkbox"/> Obsidian All-Ceramic |
| <input type="checkbox"/> PrismaTik CZ-Tuff | <input type="checkbox"/> Obsidian Fused to Non-Precious |
| | <input type="checkbox"/> Obsidian Fused to Noble |
| | <input type="checkbox"/> Obsidian Fused to WHN |

SCREW-RETAINED RESTORATIONS

- | | |
|--|---|
| <input type="checkbox"/> BruXZir Full-Strength* (w/ Ti-Base) | <input type="checkbox"/> Obsidian Fused to Noble |
| <input type="checkbox"/> NEW! BruXZir Esthetic (w/ Ti-Base) | <input type="checkbox"/> Obsidian Fused to WHN |
| <input type="checkbox"/> IPS e.max (w/ Ti-Base) | <input type="checkbox"/> PrismaTik Clinical Zirconia (w/ Ti-Base) |
| <input type="checkbox"/> Obsidian All-Ceramic (w/ Ti-Base) | <input type="checkbox"/> PrismaTik CZ-Tuff (w/ Ti-Base) |


CONTOUR AND OCCLUSION DESIGN

Embrasures: Closed* Open
Occlusion: Light* Ideal Open _____ mm Out
Contacts: Broad & Tight* Pinpoint Light

RESTORATION PONTIC DESIGN



RESTORATION METAL DESIGN



*Standard unless specified otherwise

IF NO OCCLUSAL CLEARANCE

- Call doctor Spot opposing Metal occlusion
 Metal island Make this a permanent note in my master file

Signature _____ (see reverse for limited warranty details)

License # _____ Date _____

IMPLANT WORKING TIMES

Please allow full working time for each product selected. Working times are **NOT** guaranteed and do **NOT** include weekends or holidays. Rush service available on most products but must be prescheduled (see below).

Custom Abutments and Crowns	<u>Days in Lab</u>
Custom abutment.....	8
Obsidian/BruXZir/IPS e.max/Prismatik CZ/Full-cast restoration.....	5
Obsidian/IPS e.max/Prismatik CZ/Full-cast over stock abutment.....	7
BruXZir over custom abutment.....	5
BruXZir over stock implant abutment.....	7
BruXZir, IPS e.max, or Obsidian screw-retained restoration	8
Overdentures and Fixed Dentures	
Custom impression tray	2
Bite block	5
BruXZir Full-Arch or implant denture transfer/wax setup try-in.....	5
BruXZir Full-Arch or implant denture wax setup reset.....	4
BruXZir Full-Arch or implant denture wax setup to finish.....	6
Inclusive CAD/CAM-Milled Implant Bar	7
Implant verification jig	4
Wax setup try-in, implant verification jig and custom tray	11
Final BruXZir Full-Arch Implant Prosthesis	8
BioTemps Provisionals	
BioTemps Provisionals	5
BioTemps cement over implant/cast-metal substructure/ screw-retained over implant.....	6
Custom Healing Components	
Custom healing abutment with impression coping.....	3
Custom temporary abutment with impression coping.....	5
Prosthetic stent	3



**All Restorations
Made in the USA**

*All rush cases must be prescheduled by calling **800-944-7874** before the case is shipped. Time of pickup and delivery may affect turnaround time.*

TERMS AND WARRANTY INFORMATION

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell Laboratories is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewelldental.com/policies-and-warranties.

CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS ARE COMPATIBLE WITH THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™
Certain®*

CAMLOG®
SCREW-LINE*

Dentium®
Implantium®†
SimpleLine® II†
SuperLine®†

DENTSPLY Implants
ANKYLOS® C/X*
ASTRA TECH Implant System®*

Glidewell Direct
Hahn™ Tapered
Implant System†
Inclusive® Tapered
Implant System*

HIOSSSEN®
HG System*

Keystone Dental
PrimaConnex®†

MegaGen
AnyRidge®
Implant System*

Neoss®
Neoss®
Implant System†

Nobel Biocare
Brånemark System® RP*
NobelActive®*
NobelReplace®*

Straumann®
Bone Level*
Tissue Level*

Sweden & Martina
Premium†
Shelta†

Zimmer Dental
Screw-Vent®*

*Manufactured using Inclusive® components. †Manufactured using original equipment manufacturer (OEM) components. Inclusive is a registered trademark of Prismatik Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prismatik Dentalcraft, Inc. All other trademarks are property of their respective owners. For Dentium implant systems, only titanium custom abutments are available.